

PEDIATRIC TABLES



New Design!



lb.-oz. readout



A. Model 4941

New Product!



B. Model 4944

A. Model 4941 Econo-Line™ Pediatric Table



SPECIAL FEATURES:

- New (2) polystyrene drawers and fixed shelf.
- Folkstone Gray laminate.
- 3mm high-impact ABS edges on both drawer fronts.

MODEL	LENGTH	DEPTH	HEIGHT
4941	43"	24"	36"

B. Model 4944 Pediatric Table



SPECIAL FEATURES:

- Cabinet has (2) drawers and (2) doors.
- Folkstone Gray cabinet color.
- 3mm high-impact ABS edges on doors & drawer fronts.

MODEL	LENGTH	DEPTH	HEIGHT
4944	43"	24"	36"

C. Model 4945 Pediatric Table






SPECIAL FEATURES:

- Cabinet has (2) drawers and (2) doors.
- Continuous black full-length pulls.
- Folkstone Gray or choice of any high pressure laminate solid color or woodgrain.

MODEL	LENGTH	DEPTH	HEIGHT
4945	43"	24"	36"

Common Features

PEDIATRIC TABLE WITH BUILT-IN SCALE — MODELS 4941, 4944, 4945:

- Wood has no added urea-formaldehyde for safer indoor quality.   
- PVC-Free construction.
- Digital scale with lb-oz or kg readout.
- Reads in 1 oz. increments.
- Electric and battery operation. 9V AC adapter included. Auto power shot-off. Batteries not included.
- 39" long built-in infantometer with 10" long sliding fold-away flap.
- Comfortably padded top has a protective rim on (3) sides.
- Built-in Paper Dispenser and Cutter. Holds 14" wide table paper.

COLORS: (FOLD OUT FRONT COVER)

Upholstery – 5 PVC-Free Urethane colors.



C. Model 4945