

## A. Model 5556 Mobile Therapy Station



### FEATURES:

- Durable laminate finish.
- Top laminate surface is 36" W x 18"D.
- 3mm PVC edging on all exposed edges.
- (5) High Impact polystyrene 18" W drawers.
- 18"W Oak Laminate Door with Lock.
- Grab-bar and durable 4" locking casters for mobile use.
- Hooks for (10) dumbbells.
- Almond coated pegboard with (10) 2" safety hooks for Hugger® Weights.

*New Product!*



A. Model 5556  
Front view



Model 5556  
Reverse view

MODEL #	LENGTH	DEPTH	HEIGHT
5556	44"	24"	42"

## Model 5556-100 Accessorized Rack



### CONSISTS OF

- (1) 5556 Rack, (2) 5503 Sets of (5) Dumbbells,
- (1) 5541 Set of (5) Weight Bars, (2) 5588 Sets of (5) Hugger® Weights.

## B. Model 5561-100 Pulley/Mirror Rack



- Built-in #2629 Duplex Pulleys (see page 55).
- Full view glass mirror with ANSI safety backing.
- (2) Dumbbell Racks including (20) Dumbbells.
- Durable black laminate finish on all surfaces.
- Base plate has bumper T-molding and (4) 4" dia. swiveling locking casters.

MODEL	WIDTH	DEPTH	HEIGHT
5561-100	24"	28"	73"



B. Model 5561-100  
Front view



Model 5561-100  
Reverse view

## C. Model 5594 Exercise/Rack Station



- Built-in #2629 Duplex Pulleys. (See pg. 55)
- Vertical Dispenser for (5) REP Bands™
- Storage Bin for (6) "Soft-Ball" Weights.
- Storage for Weight Bars & (3) Plastic Tubs.
- Pegboard for Hugger® Weights with (11) hooks.
- (20) Dumbbell hooks.
- (4) 4" diameter swiveling locking casters.

MODEL	WIDTH	DEPTH	HEIGHT
5594	32"	37"	74"

## Model 5594-100 Accessorized Rack

MODEL	CONSISTS OF
5594-100	(1) 5594 Rack with accessories as follows: (1) 5581 Set of (16) Hugger® Weights, (2) 5505 Sets of (10) Dumbbells ea., (1) 5531 Set of (5) REP Bands™, (1) 5541 Set of (5) Weight Bars, (1) 5545 Set of (6) "Soft-Ball" Weights

*Weights, dumbbells and other accessories must be purchased separately for Models 5556 & 5594.*



C. Model 5594  
Front view



Model 5594  
Reverse view