

**B. Optional Model 2525  
Finger Ladder**



**A. Model 6690**



View with easel side up

**A. Model 6690  
Cubex™ Therapy System  
on Wheels**



**FEATURES:**

- Adjustable height 27" to 37" foldaway 90° bolt-board work surface on right side. Includes 24 assorted bolts and (2) screwdrivers, (2) wrenches, (1) Allen wrench.
- Adjustable height 27" to 37" foldaway work surface with pegboard pad on left side. Elevates to 60° or 90°.
- Folkstone Gray laminate.
- Durable dowel and glue construction.
- 4" twin-wheel locking casters.
- Includes (100) plastic stackable pegs.
- Slide-out laminate work surface in rear.
- Reversible Activity Board in top has pegboard pad on one side and smooth easel on the other side. Board rotates 360° and locks in flat or at 60° angle.
- Two door cabinet with lock and (4) tubs.

MODEL	LENGTH	DEPTH	HEIGHT
6690	38"	24"	38"

**B. Model 2525  
Table Mounted Finger Ladder**



**FEATURES:**

- (15) peg holes in each of (3) sides.
- (50) plastic stackable pegs and C-Clamp.

MODEL	LENGTH	DEPTH	HEIGHT
2525	2½"	2½"	34"

**C. Model 6695  
Mobile Cabinet for  
Splinting and Supplies**



**FEATURES:**

- Grab-bar handle for mobility on left side.
- Door cabinet on right with lock.
- High impact polystyrene drawer.
- (3) vertical cubby sections for splinting material behind hinged door.
- Folkstone Gray laminate.
- Durable dowel and glue construction.
- 4" twin-wheel locking casters.
- Splinting pan is not included.

MODEL	LENGTH	DEPTH	HEIGHT
6695	42"	24"	31"



**C. Model 6695  
(splinting pan  
not included)**