

PORTABLE TREATMENT / SIDELINE TABLE **PROTEAM** by HAUSMANN™



600
LBS

Model 7650



Face Cradle with Face Cushion



Turf Pad

Model 7650

Portable Treatment/Sideline Table

- 2 ½" Thick multi-layered semi-firm foam.
- Comfortable PVC-Free urethane upholstery.
- Lightweight aluminum frame weighs only 30 lbs.
- Push button height adjustment.
- Includes adjustable face cradle with ultra soft face cushion and carry bag with shoulder sling.
- Spacious 73" long x 28" wide top.
- Height adjustable from 24" to 32".
- Removable 2" Turf Pad (Set of 4) included as a courtesy for outdoor use.
- Colors: 751 Black Standard (752 Blue available as special order).



Carry Bag

MODEL	LENGTH	WIDTH	HEIGHT
7650	73"	28"	24" - 32"