

COMMON FEATURES

SPLIT LEG TABLES - MODELS A9063, A9088, A9718:

- Extremely versatile and functional split leg design.
- Each padded leg rest can be adjusted and locked into desired angle up to 45° with spring-loaded lock-in lever.
- Comfortable 2" thick high-density foam upholstered tops.

COLORS: (SEE INSIDE FRONT COVER)

Upholstery - 10 Pro-Form™ vinyl colors



A. MODEL A9063

Shown with option:
69 - Gas Spring Backrest

**A. MODEL A9063
SPLIT LEG TABLE**



- Sliding doors on one side. Solid finished back.
- Natural Oak, Wild Cherry or Folkstone Gray laminate.

MODEL	LENGTH	WIDTH	HEIGHT
A9063	78"	30"	31"

**B. MODEL A9088
6-LEG SPLIT LEG TABLE**



- Six legs and full plywood subtop for added strength.
- Laminate shelf with (2) long stretchers.
- High-pressure laminate legs, apron and H-Brace.
- New enhanced leg hardware system.
- Legs have protective black vinyl scuff caps.
- Natural Oak or Wild Cherry laminate.

MODEL	LENGTH	WIDTH	HEIGHT
A9088	78"	30"	32"



B. MODEL A9088

C. MODEL A9718



**C. MODEL A9718
ELECTRIC SPLIT LEG TREATMENT TABLE**




- Durable motorized operation with "hands-free" foot control.
- Electro-hydraulic cylinders are built into (4) black powder coated steel legs with levelers.
- Fully enclosed oak laminate cabinet with sliding laminate doors.

MODEL	LENGTH	WIDTH	HEIGHT
A9718	78"	30"	27" - 37"

OPTIONS

MODELS A9063, A9088, A9718

*69 - Gas Spring Backrest 0° to 75° 

*99-L - Left hand table (A9063 and A9718.)

*FACTORY INSTALLED