

**A. MODEL A9021
RECOVERY TABLE**



- Foot end of table has a gas-spring elevating leg rest (up to 75°) and a full-extension modality drawer.
- Middle section of table has a side pullout drawer with (4) drying rods.
- Head section of table has (2) hinged doors with interior adjustable shelf offering concealed storage.
- Finished laminate back.
- Table can be used as a traditional treatment table with elevating backrest by reversing the athlete.
- Grommet holes in foot section of table and in black laminate base allow for easy cord management.
- 2" thick high-density foam upholstered tops.
- 3mm high-impact PVC edges on all door and drawer fronts.

COLORS: (SEE INSIDE FRONT COVER)

Laminate - Natural Oak and Wild Cherry

Upholstery - 10 Pro-Form™ vinyl colors

MODEL	LENGTH	WIDTH	HEIGHT
A9021	78"	30"	31"

A. MODEL A9021

Shown as right hand table (standard configuration)
Shown with option: 32 - Pillow



OPTIONS

MODELS A9021, A9075, A9078

- *A945 - Drawer Lock
- *A946 - (1) Lock/Bolt for Double Doors
- *99L - Left Hand Table

*FACTORY INSTALLED



Scan QR code for more information

B. PRO-OPTIONS TABLE



- Designed to allow you to customize to meet your specific needs... it is as easy as 1-2-3.
- Choose your cabinet style - open cabinets on both sides of table or cabinets on one side of table with solid back.
- Add optional elevating split leg feature and/or gas-spring backrest.
- Customize each compartment with any of five possible options (some restrictions).
- 2" thick high-density foam upholstered tops.
- 3mm high-impact PVC edges on all doors and drawer fronts.

COLORS: (SEE INSIDE FRONT COVER)

Laminate - Natural Oak and Wild Cherry

Upholstery - 10 Pro-Form™ vinyl colors

MODEL	LENGTH	WIDTH	HEIGHT
A9075	78"	30"	31"
A9078	78"	30"	31"

B. MODEL A9075

Shown with options:
A962 - (2) Doors, A963 - Pass thru drawer
A966 - (3) Drawers, 69 - Gas Spring Backrest, A961 - Split Leg



- STEP 1 - Choose your cabinet style**
(standard configuration is flat top):
- A9075 - Open Cabinet Both Sides†
 - A9078 - Cabinet with Solid Back

STEP 2 - Choose your top options:



69 - Gas Spring Backrest, 0° to 75°



A961 - Elevating Split Legs

STEP 3 - Customize your compartments with any of the following options
(standard compartment configuration is open with (1) adjustable shelf)



- A962 - Double Doors (includes (1) adjustable shelf)
- A963 - Pass-Thru Modality Drawer (A9075 only - includes 2 grommets for electrical cord management)
- A964 - Modality Drawer (A9078 only - includes 2 grommets for electrical cord management)
- A965 - 2-Drawers†
- A966 - 3 Drawers†

† NOTES (for A9075 Open Cabinet Table)

- Head section of table has offset support (not open pass-thru compartment - only option available is A962 Double Doors).
- If A965 or A966 Drawer option is selected, opposite side will be solid panel.