


COMMON FEATURES

PROTEAM MODULAR TAPING STATIONS:

- Upholstered Taping Modules are 36"H or 42"H. Storage Modules with laminate tops are 4" lower.
- You can order individual taping or storage modules. See page 6 or go to www.proteamtapes.com.
- Durable dowel and glue laminated wood construction for long term structural integrity.
- Self-closing laminate drawer with nylon rollers.
- Doors have concealed, adjustable, self-closing 110° Euro-hinges.
- All doors and drawer fronts have 3mm high impact edging.
- Cabinet backs are finished.
- Taping station end units have side open storage with (1) adjustable shelf.
- Comfortable 2" thick high-density foam upholstered tops are removable without tools.



COLORS: (SEE INSIDE FRONT COVER)

Laminate* - Natural Oak, Wild Cherry, or Folkstone Gray

Upholstery* - 10 Pro-Form™ vinyl colors

*Custom laminate and upholstery available

A. MODEL A9521

with options



A. MODEL A9521

2-SEAT TAPING STATION WITH OPTIONS

SHOWN L TO R:

- A918 Full Extension Pullout
- A926 Fixed 75° Upholstered Back
- A923 Set of (3) Drawers
- A924 Waste Hamper
- A945 Top Drawer Lock
- A944 Full Height Doors
- A946 Lock/Bolt
- A929 Side Pull-Out Stand

MODEL	LENGTH	DEPTH	HEIGHT
A9521	6'	42"	36"

B. MODEL A9400

SINGLE SEAT TAPING STATION

- Full extension pullout for tape storage opens to the side.
- Two 8" self-closing laminate drawers with nylon rollers for extra deep storage to the side for easy access to supplies.
- Large side storage compartment with shelf includes door with concealed, adjustable, self-closing 110° Euro-hinges.
- Left hand configuration is standard.
- Taping station can be configured left or right. Please specify for right hand configuration

MODEL	WIDTH	DEPTH	HEIGHT
A9400	28"	42"	36"
A9401	28"	42"	42"



B. MODEL A9400

Shown as left hand unit (standard configuration)

Option A926
Fixed 75°
Upholstered Back



OPTIONS

MODELS A9400, A9401

- *A926 - Fixed 75° Upholstered Back
- *A945 - Door/Drawer Lock
- *99-R - Pullout, drawers, and door are on left when facing taping station

*FACTORY INSTALLED