

Common Features

**PARALLEL BARS —
MODELS 1318, 1319, S-320, S-321:**

- 1½" diameter one-piece stainless steel handrails for easy use and durability.
- Heavy gauge steel base plates, uprights and fittings. Black powder coated chip and wear resistant finish.
- Each upright telescopes up in 1½" increments and locks into (10) height positions with fail-safe ball-tip locking pin.

A. Models 1318, 1319 Folding Parallel Bars



SPECIAL FEATURES:

- Designed to easily fold to save space.
- Unit is only 10" wide when folded.
- Height adjustable from 28" to 41".
- Fixed width dimension of 25".

MODEL	LENGTH	WIDTH	HEIGHT
1319	7'	25"	28"- 41"
1318 (not shown)	9'	25"	28"- 41"

B. Models S-320, S-321 Econo Platform Mounted Parallel Bars



SPECIAL FEATURES:

- Folkstone Gray laminate platform with tapered ends and anti-slip rubber treads.
- Height and width adjustable bars.

MODEL	LENGTH	WIDTH	HEIGHT
S-320	10'	15"- 28"	28"- 41"
S-321 (not shown)	7'	15"- 28"	28"- 41"

*platform is 36" wide



A. Model 1319

Model 1319
(folded)



B. Model S-320