



## Advantages

- EMR/EHR Connectivity Featuring USB Port**  
Scale has the ability to interface with a computer, monitor, printer or other electronic device via USB
- Hospital Grade Construction**  
Durable steel structure with large sturdy stretcher enables patients to be weighed accurately and comfortably
- Flexible Functionality**  
4 wheels for ease of movement, 2 locking wheels for safety, 70" x 30" stretcher allows for efficient operation
- Value-Added Features**  
Operating functions: BMI, Tare, Zero Out, Hold, Live Weight



2 Locking wheels  
for safety

Model	Master Pack	Shipping Dims. (L x W x H) / Cube (3 cartons)	Weight (LB) total: 142 lb	UPC Code
2000KL	Stretcher	71.75" x 3.5" x 2.5" / 0.36	10-1/2 lb	892076002466
	Lifter	51.5" x 25.75" x 13" / 9.98	86-1/2 lb	
	Electronics			



**Specifications**

- Capacity: 400 lb / 181 kg
- Resolution: 0.2 lb / 0.1 kg
- Connectivity via USB
- 3/4" LCD display
- Functions: LB / KG Conversion, LB / KG lock out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Reweigh, Auto Zero, Auto Off
- Stretcher: 70-3/8" (w) x 30-3/8" (d)
- (4) Wheels for ease of movement, 2 locking wheels for safety
- 120V adapter included, 6 AA batteries (not included)
- 2 Year Limited Warranty



2000KL: Reduced or minimum foam in packaging / recyclable cardboard packaging.