

## THE 221 HR

# TAA COMPLIANT PORTABLE STADIOMETER WITH CARRYING CASE

The number-one-selling medical scale company in North America provides the optimum in quality, reliability, and ease of use with the 221 HR. This TAA Compliant stadiometer, designed to meet the needs of doctor's offices, schools, home health and traveling clinicians, features an easy to assemble / disassemble structure and a large sturdy anti-slip base. The long headpiece adds comfort for patients of all sizes and the large view finder provides an easy-to-read accurate measurement. The 221 HR has a wide measuring range of 17.7 cm - 210 cm (7" - 82.6"). Industry leading 1/8" platform height, requires a minimal step up, offering the patient a more comfortable and stable weighing experience. This light weight stadiometer includes a plastic carrying case and gives clinicians a cost effective option to measure all types of patients. It can

be assembled and disassembled in seconds and is free-standing as wall fastening is not required.



### Large View Finder

The large view finder moves up and down providing an easy-to-read accurate measurement.

### Carrying Case Included

The 221 HR includes a polycarbonate carrying case which makes the stadiometer easy to move from location to location.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers to our stadiometer, it's not surprising that the 221 HR is accurate, easy to use, portable and convenient.

**150,000+**  
STADIOMETERS  
INSTALLED  
NATIONWIDE

**120+** LOCAL  
FIELD SPECIALISTS


**1**  
CARRYING CASE  
INCLUDED

**#1** TRUSTED  
BRAND

**100** YEARS  
OF SERVICE

**Health o meter**<sup>®</sup>  
Professional

Health o meter<sup>®</sup> Professional 221 HR



weigh**easier**<sup>®</sup>

# TECHNICAL SPECIFICATIONS

### Measuring Range

7" – 82.6" / 17.7 cm – 210 cm

### Graduation

1/16" / 1 mm

### Platform Size (w x d x h)

13.75" x 14.25" x 0.125"  
(406 mm x 365 mm x 3mm)

### Product Footprint (w x d x h)

13.75" x 23" x 89"  
(349 mm x 585 mm x 2261 mm)

### Product Weight (Stadiometer)

7.9 lb / 3.6 kg

### Product Weight

(Stadiometer in case)

10.2 lb / 4.6 kg

### Carrying Case

Plastic carrying case included

### Carrying Case Size (w x d x h)

25" x 5" x 15"  
(635mm x 127mm x 381 mm)

### Warranty


2-year limited manufacturer warranty

# ORDERING INFORMATION

MODEL	DESCRIPTION
221HR	Portable stadiometer with carrying case




Carrying Case Included



Anti-Slip Base

Industry Leading 1/8" Platform Height



To learn more, call **1-800-815-6615**  
or visit [hom scales.com/221HR](http://hom scales.com/221HR)

**Health o meter**<sup>®</sup>  
Professional

