



## Advantages

- ✓ **EMR/EHR Connectivity Featuring USB Port**  
Scale has the ability to interface with a computer, monitor, printer or other electronic device via USB
- ✓ **Hospital Grade Construction**  
Durable steel structure with non-slip platform enables all size patients and wheelchairs to be weighed accurately and comfortably
- ✓ **Value-Added Features**  
Operating functions: BMI, Tare, Zero Out, Hold, Display Adjustment, Live Weight
- ✓ **Dual Ramps and Remote Display**  
Non-slip ramps provide easy on and off access to weighing platform and display can be mounted on wall or set on table allowing for convenient placement



**2650KL** Platform: 35-3/4" (w) x 32-1/4" (d)

**2700KL** Extra Large Platform: 43" (w) x 42" (d)

Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2650KL	1			892076002442
2700KL	1			892076002633

2650KL

2700KL



**Specifications**

- Capacity: 1000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Connectivity via USB
- 3/4" LCD display, can be wall mounted or placed on flat surface
- Functions: LB / KG Conversion, LB / KG lockout, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Reweigh, Auto Zero, Auto Off
- Platform: 35-3/4" (w) x 32-1/4" (d) (2650KL)
- Platform: 43" (w) x 42" (d) (2700KL)
- Ramp (dual) Size: 10" x 32-1/4" (2650KL)
- Ramp (dual) Size: 12-3/4" x 42" (2700KL)
- (2) Wheels for ease of movement
- 120V adapter included, 6 AA batteries (not included)
- 2 Year Limited Warranty



2650KL & 2700KL: Reduced or minimum foam in packaging / recyclable cardboard packaging.