



Advantages

- Stainless Steel Base and Weighing Surface**
The stainless steel design is easy to clean and ideal for sterile environments.
- Rechargeable Battery Power Source**
Scale is powered by a 6V rechargeable battery with AC adapter/battery charger included with the scale.
- High Resolution Measurements**
Scale provides highly precise and accurate results with refined 0.05 oz / 1 g resolution.
- Optional Tray**
9 ¾" x 7 ¾" tray is useful for weighing items that require a flat surface. Stainless steel tray is removable, can be autoclaved, and is available for purchase in packages of 3.
- Convenient Weighing Pan**
The 2 ½" deep and 10" wide stainless steel weighing pan effectively retains fluids and can be removed for sterilization. Additional pans can be purchased in packages of 3.
- Tare Functionality**
Container weight can be easily deducted to provide accurate net weight of the contents.
- Backlight Display**
Easy to read from any angle and light condition.



Optional Weighing Tray Available
Order # 3400TRAY



Additional Weighing Pans Available Order # 3400PAN

| Model | Master Pack | Shipping Dims. (L x W x H) / Cube | Weight (LB) | UPC Code |
|--------|-------------|-----------------------------------|-------------|--------------|
| 3400KL | 1 | 14" x 14" x 11.5" / 1.3 | 13.7 lb | 854896003788 |



Specifications

- Capacity: 176 oz / 5000 g
- Resolution: 0.05 oz / 1g
- Weighing Surface Dimensions (pan): 10" (w) x 9-7/8" (d) x 2 1/2" (h) / 254mm (w) x 249mm (d) x 61mm (h)
- Weighing Surface Dimensions (optional tray): 7-3/4" (w) x 9-3/4" (d) x 5/8" (h) / 95.25mm (w) x 248mm (d) x 15.875mm (h)
- Scale Base Dimensions (including handles): 11 1/2" (w) x 9 7/8" (d) x 5 1/4" (h) / 292mm (w) x 250mm (d) x 133mm (h)
- Display: 5/8" (15.9mm) LCD Display
- Power Source: 100 – 240 VAC AC adapter/battery charger / 6 volt rechargeable battery (included)
- Functions: Zero Out/Tare, Auto Zero, Auto Off, g/oz Conversion, Low Battery Indicator, Backlight Display
- 2 Year Limited Warranty

Accessories Available

| | |
|-----------------|--|
| 3400PAN | Optional 3 pack of stainless steel pans |
| 3400TRAY | Optional 3 pack of stainless steel trays |