

THE 500KL

PRIMARY CARE DIGITAL SCALE

The number-one-selling physician scale in North America, the all-digital Health o meter® 500KL primary care scale provides the optimum in quality, reliability, and ease of use. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale features a built-in height rod in the scale column, a large non-skid platform, and wheels that enable portability. The scale gives providers and patients an eye-level weight indicator for patients up to 550 pounds, with a user-friendly digital interface that includes BMI, hold/release, zero, tare, reweigh, recall, KG/LB, clear, and enter functions. The 500KL offers the best in its class extra-wide 8" wheelbase making transport easy. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The first healthcare scale to provide easy integration and seamless communication with today's EMR applications, the 500KL includes the option to add wireless connectivity to interface with a computer, EMR software, or other electronic device.

With its low platform height, the 500KL is easy to set up and offers patients a more comfortable and stable weighing experience.

500HB

"live" wrap-around handlebar

The 500HB is a unique handlebar accessory for the 500KL and is easy to attach to any new scale using the pre-drilled mounting holes and can also be attached to any existing scale using the easy-to-use drill guides. It provides the stability that many patients need without affecting the accuracy of the weight reading.

THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 500KL is the most widely used scale in physician offices nationwide.

150,000+
INSTALLED
NATIONWIDE

120+ LOCAL
FIELD SPECIALISTS

8" WIDE WHEELBASE
BEST IN ITS CLASS


#1 TRUSTED
PHYSICIAN SCALE

100 YEARS
OF SERVICE

Health o meter®
Professional

Health o meter® Professional 500KL

weigh**easier**®



TECHNICAL SPECIFICATIONS

Capacity

550 lb / 250 kg

Resolution

0.2 lb / 0.1 kg

EMR connectivity

Via optional Pelstar® wireless technology

Display

1" / 25 mm, LCD display

Batteries

6 AA batteries included, optional adapter (order# ADPT31)

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare, Pre-tare
Hold/Release
Auto Zero
Auto Off
Reweigh
Recall
Parent/Child weighing function
Everlock® - permanent unit lock

Platform Size (w x d x h)

13 3/4" x 16 1/2" x 2 3/8"
(349 mm x 419 mm x 60 mm)
Dent and rust proof platform

Product Footprint

14 3/8" x 21 1/8" x 53 5/8"
(359 mm x 537 mm x 1362 mm)

Product Weight

21 lb / 9 kg

Height Rod

30"-84 1/4" / 76 cm-214 cm
Graduation: 1/8" / 1 mm

Portability

2 wheels, wide 8" wheelbase for ease of movement

Warranty

2-year limited manufacturer warranty



ORDERING INFORMATION

MODEL	DESCRIPTION
500KL	Eye-level digital scale, LB & KG
500KG	Eye-level digital scale, KG only
500KLHB	500KL with handlebar accessory
500KLAD	500KL with adapter ADPT31
500KGAD	500KG with adapter ADPT31
500KL-BT	500KL with built-in Pelstar® wireless technology
500KG-BT	500KG with built-in Pelstar® wireless technology
500HB	Optional handlebar accessory for 500 series scales
SS-500KL	ScaleSurance® 2 year extended warranty for 500 series*
ADPT31	Optional adapter for 500KL & 500KG
ADPT30	Optional adapter for 500KL-BT & 500KG-BT
PLUGSET30	Plug adapters for use outside North America

To learn more, call **1-800-815-6615**
or visit homscapes.com/500KL

Health o meter®
Professional



*Not available in all countries.