

## **Advantages**

**▼** EMR/EHR Connectivity via USB

Scale has the ability to interface with a computer, monitor, printer or other electronic device via USB for fast, easy and accurate data transfer

Sanitary Plastic Weighing Surface

The easy-to-clean plastic tray is ideal for sterile environments

**Built-In Measuring Tape** 

Easy to read measuring tape allows user to measure both weight and height simultaneously **▼ Display Head 17" Above Scale Base** 

Display head is positioned 17" above the scale base to allow for an unobstructed view of the weight measurement

**✓ Value-Added Features** 

Operating functions (LB/KG Lockout, Tare, Zero Out, Hold/Release, Auto Zero, Auto Off) included at no additional cost





Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
522KL	1	17.13" x 26.25" x 14.25" / 3.71	23 lb	892076002541



## **Specifications**

- Capacity: 50 lb / 23 kg
- Resolution: 0 20 lb 0.2 oz / 20 50 lb 0.5 oz 0 9 kg / 5g; 9 23 kg / 10 g
- EMR connectivity via USB
- 1" LCD display
- Functions: LB / KG conversion, LB /KG lock out, Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Tray Size: 24 1/16" (w) x 14 1/8" (d) x 3" (h)
- Measuring Tape: 0" 23" / 0 cm 58 cm
- 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- 2 Year Limited Warranty

