




Remote Display can be mounted on a wall or placed on a table

Advantages

- EMR/EHR Connectivity via USB**
Scale has the ability to interface with a computer, monitor, printer or other electronic device via USB for fast, easy and accurate data transfer
- Body Mass Index (BMI) Capability**
Scale functions include the ability to automatically calculate the patient's Body Mass Index (BMI)
- Remote Display**
The remote display provides convenience with the ability to be mounted on a wall or placed on a table
- Value-Added Features**
Operating functions (BMI, LB/KG Lockout, Hold/Release, Auto Zero, Tare, Zero Out, Auto Off) included at no additional cost
- 600 lb Capacity Floor Scale**
The 752KL is the market's highest capacity floor scale, providing more versatile use
- AC Adapter Included**
Scale is powered by a 120V AC adaptor (included) or 6 C-cell batteries (not included) which allows for the convenience of dual power sources



PORTROD
Optional Wall Mounted
Height Rod Available



Remote Display can be mounted on a wall or placed on a table

Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
752KL	1	24.4" x 16.14" x 4.53" / 1.03	18 lb	892076002138

Specifications



- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- EMR connectivity via USB
- RJ45 Serial Port for printing
- Easy to read 1" LCD display
- Remote Display can be placed on table, mounted on a wall or held by the operator
- Functions: LB / KG conversion, LB /KG lock out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Platform Size: 14-3/4" (w) x 14-11/16" (d)
- 120V adapter included, 6 C-cell batteries not included
- Optional Wall Mounted Height Rod available (item #PORTROD)
- 2 Year Limited Warranty



752KL has no foam in packing and recyclable cardboard packaging