



Advantages

- Choice of Units of Measurement**
Digital display provides height measurement in centimeters, feet, or inches.
- Portability**
Lightweight and battery powered, the height rod can be moved anywhere it is needed.
- Efficient Design**
Measuring infant length is fast and easy with clear digital measurements, large head and foot paddles, and convenient Hold function.

Specifications

- Measurement Range: 13 ¾" – 31 ½" / 35 -80 cm
- Resolution: 1/16" / 0.1 cm
- Physical Dimensions:
 - Depth: 11 ½" / 290 mm
 - Width (when closed): 24 ¾" / 630 mm
 - Height: 2 ¾" / 70 mm
- Power Source: 2 AAA batteries (included)
- Display: LCD, 5/8" / 16 mm digit height
- Functions: Auto Off, Hold
- 2 year limited warranty



The eBABYHR attaches to the 550KL for a comprehensive weighing and measuring device

Model	Description
SS-eBABYHR	ScaleSurance 2 year extended warranty*