

ELEVATE
The Health o meter[®] EMRscale



Advantages

- ✓ **EMR/EHR Connectivity via USB**
Scale has the ability to interface with a computer, monitor, printer or other electronic device via USB for fast, easy and accurate data transfer
- ✓ **350° Swivel Display**
Eye-level rotating display can be read from almost any angle and allows for patient privacy and ease of use
- ✓ **Digital Height Rod**
Fast and easy-to-use digital height rod provides option to capture height measurement electronically and accurately
- ✓ **Value-Added Features**
Operating functions (BMI, LB/KG Lockout, Hold/Release, Auto Zero, Tare, Zero Out, Auto Off) included at no additional cost
- ✓ **Body Mass Index (BMI) Capability**
Scale functions include the ability to automatically calculate the patient's Body Mass Index (BMI)
- ✓ **EMR/EHR Compliant**
The only scale available that meets the new government healthcare regulations regarding BMI and EMR/EHR compliance

Heavy Duty Eye Level Digital Scale with Integral Digital Height Rod

ELEVATE EMRscale™

Health o meter®
P r o f e s s i o n a l

ELEVATE
The Health o meter® EMRscale



Model : ELEVATE-C

To be used with the Welch Allyn
Connex® Vital Signs Monitors

Model : ELEVATE-L

To be used with the Welch Allyn
Spot Vital Signs® LXi

Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
ELEVATE-C	1	52" x 26" x 6.3" / 4.5	38	854896003085
ELEVATE-L	1	52" x 26" x 6.3" / 4.5	38	854896003092

Specifications



- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- EMR connectivity via USB
- 1" LCD display, 35° left / right swivel (350° with height rod fully retracted) for ease of viewing
- Functions: LB / KG conversion, LB /KG lock out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- Platform Size: 2-1/2" (h) 14-3/4" (w) x 14-11/16" (d)
- Digital Height Rod: 24" - 88" / 61 cm - 223 cm
- (2) Wheels for ease of movement
- 120V adapter included, 6 C-cell batteries not included
- 2 Year Limited Warranty



ELEVATE has no foam in packing and recyclable cardboard packaging

Health o meter®
Professional
www.homscales.com

We Always Accept:
VISA, MC, & AMEX



Health o meter® Professional Scales
11800 South Austin Avenue
Alsip, IL 60803 USA
Phone: 800-815-6615
Fax: 800-654-7330

**BACKLOG
FREE ZONE**

We will do everything possible to make it easy for our customers to do business with us.