



Staph

E. coli

lifeSign

INFECTIOUS DISEASE PRODUCTS



Strep

MRSA

Strep Grouping Kit (Enzyme)

A latex streptococcal grouping for rapid ID of Beta-hemolytic streptococcal groups A,B,C,D,F & G from culture media



Simple Procedure

Efficient enzyme extraction procedure
Easy-to-use latex agglutination format
Easy-to-read blue latex on a white card

Convenient

Color-coded components
Individual reagents available separately
Space-saving work station
Results in 12 minutes

Superior Support

Knowledgeable sales staff
Responsive technical service team

ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #	MFG #
Streptococcal Grouping Kit	50	50040	(DR0585A)
Streptococcal Grouping Reagent A	50	50041	(DR0586G)
Streptococcal Grouping Reagent B	50	50042	(DR0587G)
Streptococcal Grouping Reagent C	50	50043	(DR0588G)
Streptococcal Grouping Reagent D	50	50044	(DR0589G)
Streptococcal Grouping Reagent F	50	50045	(DR0590G)
Streptococcal Grouping Reagent G	50	50046	(DR0591G)
Extraction Enzyme	50	50047	(DR0593G)
Polyvalent Positive Control	50	50048	(DR0592G)
Disposable Reaction Cards (White)	50	50049	(DR0500G)

Strep Plus Kit (Acid)

A latex streptococcal grouping for rapid ID of Beta-hemolytic streptococcal groups A,B,C,F & G from culture media

Simple Procedure

Room temperature extraction (Nitrous Acid)
Standardized latex agglutination procedure

Convenient

Ideal for "must have results"
Color coded reagents and space-saving workstation
All components available separately*
Rapid results in 2 minutes

Superior Support

Knowledgeable sales staff
Responsive technical service team



*See Grouping Reagents A-G, Control and Cards above

ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #	MFG #
Strep Plus Kit	50	50030	(DR0575M)
Strep Plus Extraction Reagent	50	50037	(DR0574M)

LifeSign PYR

Rapid colormetric test for PYRase activity in strep from culture media

Convenient

Easy to use
Large reaction area
Simple instructions printed on the test card
Rapid results in 7 minutes
Area for patient name/ID

Economical

Less expensive than "disk tests" or bile esculin and salt

More Options

Complete kit
QC or additional samples can be run concurrently



ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
LifeSign PYR	50	50060

LifeSign Staph Latex

Rapid Staph aureus test from culture media



Simple Procedure

Detects clumping factor and protein A
Rapid results in 2 minutes

Accurate:

97.6% sensitivity
98.1% specificity

Economical

Low cost test when compared to large equipment

ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
LifeSign Staph Latex	50	50051

AnaeroGen Gas Pack 2.5 L

An anaerobic atmosphere generation system

Simple Procedure

Quick and convenient system which generates an anaerobic atmosphere

Convenient

No catalyst or water is required and no hydrogen gas produced



ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
Anaerogen Gas Pack 2.5L	10	AN025

E. Coli 0157 Latex Test

A highly sensitive latex agglutination test for the identification of E. coli 0157 from culture media



Easy-to-read

Unique blue latex reagent for clearer interpretation of results

Accurate

Reliable results showing 100% sensitivity and 99% specificity

Convenient

Complete kit packaged in a space saving workstation with color coded bottle caps to easily distinguish between positive and negative controls and test reagents

Rapid results in 2 minutes

Safe

White disposable reaction cards for quick and easy disposal

ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #	MFG #
E. Coli 0157 Test	100	50052	(DR0620A)

PBP2' M-R Staph Test Kit

A latex agglutination test for the detection of MRSA from culture media

Reliable

An aid in the identification of methicillin resistance Staph aureus

Rapid

Results in 13 minutes

Convenient

Easy to use procedure
Performed on disposable reaction slides



ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
PBP2' Latex Agglutination Test	20	DR0900A

ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #	MFG #
Streptococcal Grouping Kit	50	50040	(DR0585A)
Streptococcal Grouping Reagent A	50	50041	(DR0586G)
Streptococcal Grouping Reagent B	50	50042	(DR0587G)
Streptococcal Grouping Reagent C	50	50043	(DR0588G)
Streptococcal Grouping Reagent D	50	50044	(DR0589G)
Streptococcal Grouping Reagent F	50	50045	(DR0590G)
Streptococcal Grouping Reagent G	50	50046	(DR0591G)
Extraction Enzyme	50	50047	(DR0593G)
Polyvalent Positive Control	50	50048	(DR0592G)
Disposable Reaction Cards (White)	50	50049	(DR0500G)
Strep Plus Kit	50	50030	(DR0575M)
Strep Plus Extraction Reagent	50	50037	(DR0574M)
LifeSign PYR	50	50060	
LifeSign Staph Latex	50	50051	
Anaerogen Gas Pack 2.5L	10	AN025	
E. Coli 0157 Test	100	50052	(DR0620A)
PBP2' Latex Agglutination Test	20	DR0900A	

For more information or to place an order, contact:

The logo for LifeSign features a blue curved line above the word "lifeSign". "life" is in a light blue color and "Sign" is in a green color.

71 Veronica Avenue • Somerset • New Jersey • 08873

800.526.2125 • 732.246.3366 • Fax: 732.246.0570

www.lifesignedmed.com