


3M™ Steri-Drape™  
Arthroscopy Drape with Pouch – 1194 and 1094



## Two Pouch Choices

### One-Piece Convenience


3M Steri-Drape Arthroscopy Drapes offer you the choice of two unique fluid collection pouches. Paired with 3MFabric, it helps reduce the risk of exposure to blood and body fluids for patients and healthcare workers.

- One-piece draping with bended or extended leg
- Tube and cord organizer keeps operative site clear

## 3M™ Steri-Drape™ Arthroscopy Drape with Pouch – 1194 and 1094


### 3M™ Steri-Drape™ Arthroscopy Drape with Pouch 1194

- This drape has a wide, rectangular fluid collection pouch that adapts to a bended or extended leg and allows full manipulation of the knee during surgery.
- Two elastomeric apertures isolate the operative site at the knee and ankle
- Low profile pouch stays out of surgeon's way
- Exit ports on both sides of pouch provide better fluid control
- Malleable metal wire keeps pouch open for optimal fluid collection



### 3M™ Steri-Drape™ Arthroscopy Drape with Pouch 1094

- This V-shaped pouch can be adhered to the skin in any position you want, making the pouch fit your procedure
- Elastomeric aperture isolates the operative site at the knee
- Pouch can be adhered to the skin according to surgeon preference
- Exit port provides easy drainage of fluids



### 3M Fabric

- Meets Flammability Class 1 (CFR Title 16 Part 1610)
- Meets AAMI Standard PB70:2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, Selection of surgical gowns and drapes in health care facilities
- Provides AAMI level 4 barrier performance throughout entire drape
- Low linting (INDA IST 160 1-83)
- Meets AORN 2014 Recommended Practices for materials, selection, and use of surgical drapes that represent an optimal level of practice

### Ordering Information

#### 3M™ Steri-Drape™ Arthroscopy Drape with Pouch – 1194 and 1094

Cat No.	Product Description	Packaging	Minimum Order	Pricing
1194 1194NS	89" x 120" (226cm x 304cm) 115" (292cm) Width at arm extensions	1 Each, 5 Each/Case 20 Each/Case	Case Case	Case Each
1094 1094NS	89" x 125" (226cm x 320cm)	1 Each, 20 Each/Case 30 Each/Case	Case Case	Case Each

For more information about 3M™ Surgical Drapes or other 3M Health Care products:

- » Contact your local 3M representative
- » Call the 3M Health Care Customer Helpline at 800-228-3957
- » Visit [www.3M.com/infectionprevention](http://www.3M.com/infectionprevention)



#### Infection Prevention Division 3M Health Care

2510 Conway Avenue  
St. Paul, MN 55144-1000  
U.S.A.  
1-800-228-3957  
[www.3M.com/infectionprevention](http://www.3M.com/infectionprevention)

3M and Steri-Drape are trademarks of 3M.  
©3M 2014. All rights reserved.  
70-2009-5943-8