

3M™ Steri-Drape™ Craniotomy Drape with Ioban™ 2 Incise Film – 6687



Time-saving Drape

Superior Fluid Control

3M Steri-Drape Craniotomy Drape combines fluid control and incise area with antimicrobial properties. Paired with 3M Fabric, it helps reduce the risk of exposure to blood and body fluids for patients and healthcare workers.

- 3M™ Ioban™ 2 Incise Film provides continuous antimicrobial activity and a sterile operative surface at the beginning of surgery
- The fabric absorbs spills and splatters, making multiple layering unnecessary and preventing strike-through
- Large rectangular pocket helps manage fluid collection and control

Lead
The Way

3M

3M™ Steri-Drape™ Craniotomy Drape with loban™ 2 Incise Film – 6687



A combination that helps reduce the risks of contamination of the operative site.

3M loban 2 Incise Film

- Continuous broad-spectrum antimicrobial activity in the drape adhesive
- The drape is securely adhered to the skin up to the wound edge

3M Fabric

- Meets Flammability Class 1 (CFR Title 16 Part 1610)
- Meets AAMI Standard PB70:2012 *Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, Selection of surgical gowns and drapes in health care facilities*
- Provides AAMI level 4 barrier performance throughout entire drape
- Low linting (INDA IST 160 1-83)
- Meets AORN 2014 Recommended Practices for materials, selection, and use of surgical drapes that represent an optimal level of practice



Ordering Information

3M™ Steri-Drape™ Craniotomy Drape with 3M™ loban™ 2 Incise Film – 6687

Cat No.	Product Description	Packaging	Minimum Order	Pricing Unit
6687	77" x 115" (196cm x 292cm) 100" (254cm) Width at arm extensions	1 Each 10 Each/Case	Case	Case

For more information about 3M™ Surgical Drapes or other 3M Health Care products:

- » Contact your local 3M representative
- » Call the 3M Health Care Customer Helpline at 800-228-3957
- » Visit www.3M.com/infectionprevention



Infection Prevention Division 3M Health Care

2510 Conway Avenue
St. Paul, MN 55144-1000
U.S.A.
1-800-228-3957
www.3M.com/infectionprevention

3M, Steri-Drape and loban are trademarks of 3M.
©3M 2014. All rights reserved.
70-2009-7056-7