



**Littmann**<sup>®</sup>  
Stethoscopes



Introducing:  
**A new  
class of  
clarity.**



**Assess patients more confidently with the help of the new  
3M™ Littmann® CORE Digital Stethoscope.**

Introducing our most advanced stethoscope yet. Now you can get everything you love about Littmann stethoscopes—superior acoustics, comfort and quality—with up to 40x amplification (at peak frequency, versus analog mode), active noise cancellation and mobile sound wave visualization.

**It's time to listen like never before.**

# Take listening to a whole new level.

## 3M™ Littmann® CORE Digital Stethoscope

### Comfort. Confidence.

Snap tight, soft sealing eartips conform to ears for an excellent acoustic seal and comfortable fit.

### Hear more. Way more.

Up to 40x amplification\* and active noise cancellation to reduce background noise.

*\*40x amplification based on comparison of analog versus electronically assisted amplified listening with maximum volume at the peak frequency (125 Hz).*

### Hear it the way you want to.

Quickly toggle between analog and amplified listening modes—and use the tunable diaphragms on both sides of the double-sided chestpiece to hear high and low frequency sounds.

### Discover new details with digital.

Connects to Eko software, allowing you to visualize heart sound waveforms in real time. Plus, use AI to aid in the detection of heart murmurs. (Subscription required\*\*. Learn more at [ekohealth.com/ai](https://ekohealth.com/ai))

Upgrade and listen like never before at [Littmann.com/CORE](https://Littmann.com/CORE)

## Do more than ever before.

### Capture


Save and annotate 15, 30, 60 or 120-second recordings in the secure dashboard.

### Analyze

Create a library of heart sounds to monitor the progression of a disease or for education.

### Share

Securely share recordings with colleagues.



Smart device not included.

\*\*Available in US only

©2020 3M. All rights reserved. 3M, Littmann, the L Littmann logo, the L, and Cardiology IV are marks and/or registered marks of 3M. Unauthorized use prohibited. Used under license in Canada.

Eko trademarks, including CORE, are marks and/or registered marks of Eko Devices Inc.



# Listen first.

**3M Health Care**  
3M Medical Solutions Division  
2510 Conway Ave., St. Paul, MN 55144 USA  
1-800-228-3957 | [www.littmann.com](https://www.littmann.com)

70-2011-8018-2