

# 3M Skin Stapler Systems

*Innovation is at hand.*

Continual improvement through innovation. The full line of 3M Skin Stapler Systems provides the highest degrees of performance, versatility and reliability. Each is ergonomically and functionally designed to exceed your expectations of precision and effectiveness. 3M staplers are convenient, easy to use and disposable.



2

## 3M Disposable Skin Stapler Systems

Introducing the innovative difference in wound support.

**T**he 3M™ MS Disposable Skin Stapler dispenses box-style staples versus arcuate, suture-style staples used in the 3M DS skin staplers. Both the box-style and arcuate-style “floating” staples provide several performance advantages, including optimal wound closure and support.


- The rigid design holds wound edges securely by compressing the tissue between the legs of the staple.
- “Floating” staples maintain wound edge apposition without the need to fill them with tissue; allowing for appropriate edema to occur.
- Both styles offer speed of closure, elimination of cross-hatching and ease of removal for enhanced cosmetic results as well as less wound and patient trauma.

Both MS and DS Staplers are available in three configurations: 5, 15 and 25 staples.



3

The user-friendly design of 3M™ Disposable MS and DS Skin Stapler Systems dispenses “floating” staples which allow room for edema and easy removal with less trauma.



3M MS skin staplers are slightly larger than 3M DS models, yet both are designed to be compatible with 3M Multi-Shot Skin Stapler Systems.

## 3M™ Vista Disposable

# Skin Stapler System

A new *angle* on visibility.

The 3M Vista Disposable Skin Stapler dispenses box-style staples and is the first and only skin stapler engineered to provide true staple visibility from your normal vantage point. You have a better view of your work for faster closure and optimal control. Vista skin staplers are available in four configurations: 35 regular or wide and 15 regular or wide.



- **Clear Visibility** – the

angled stapler head allows a clear view of the wound area for precise staple implant.

- **Ergonomic Feel** – the Vista skin stapler fits comfortably in your hand for enhanced control and the short trigger stroke reduces fatigue.
- **Versatile Staple Release** – staples can be released from the retractable anvil in any direction.
- **Reliable** – its simple design and quality construction assure consistent performance.




# 3M Multi-Shot Disposable Skin Stapler System

Cost-effective *versatility* for both 3M™ Disposable MS and DS Skin Stapler Systems.

The 3M Multi-Shot Skin Staplers are designed to be versatile, reliable and cost-effective wound closure tools. They combine the best of both worlds – reliability and disposability. Use either stapler alone or as a cartridge inserted in the RSH-100 reusable multi-shot skin stapler handle.



- **Quick and Easy to Use** – simply press your thumb and index finger together to implant a staple in a single motion.
- **No Waste** – selecting only the number of staples required for each incision eliminates the extra expense of discarding unused staples.
- **Reliable** – uncomplicated design delivers consistent, reliable performance.
- **Simple to Remove** – All of 3M's skin staples provide for ease of removal for enhanced cosmetic results.



# 3M™ PGX Disposable Skin Stapler System

Quality + Simplicity = *Reliability*

The quality construction and design simplicity of the 3M PGX Disposable Skin Stapler meets the highest standards of reliable performance. It provides you consistent staple implants with every use. PGX skin staplers are available in four configurations: 35 regular or wide and 15 regular or wide staples.



11

- **Precision**

**Handling** – the stapler, with its slim, lightweight design, is comfortable to use and handles like a well-balanced, precision surgical instrument.

- **Reliable** – uncomplicated design delivers consistent, reliable performance.



12

## 3M Disposable Skin Staple Removers

Help bring medical  
procedures to a  
successful close.

**T**he 3M Disposable  
Skin Staple Removers  
easily remove staples  
in the same direction  
they were implanted,  
making removal  
simple with  
little pain.



*SR-1 (tweezers-style  
staple remover)*

13



*SR-3 (scissors-style  
staple remover)*

14



# Ordering Information



15

## Disposable 3M Vista skin stapler

Catalog Number	Staple Qty/Size	Items Per Box	Boxes Per Case
3995	35 Wide	6 Staplers/Box	4
3996	15 Wide	6 Staplers/Box	4
3997	35 Regular	6 Staplers/Box	4
3998	15 Regular	6 Staplers/Box	4



16

## Disposable 3M PGX skin stapler

PGX-15R	15 Regular	6 Staplers/Box	4
PGX-15W	15 Wide	6 Staplers/Box	4
PGX-35R	35 Regular	6 Staplers/Box	4
PGX-35W	35 Wide	6 Staplers/Box	4



17

## Disposable 3M cartridges for 3M MS multi-shot skin stapler


MS-5	5 Box-style	6 Staplers/Box	4
MS-15	15 Box-style	6 Staplers/Box	4
MS-25	25 Box-style	6 Staplers/Box	4



18

## Disposable 3M cartridges for 3M DS multi-shot skin stapler

DS-5	5 Arcuate-style	12 Staplers/Box	4
DS-15	15 Arcuate-style	12 Staplers/Box	4
DS-25	25 Arcuate-style	12 Staplers/Box	4



20

## Disposable 3M skin staple removers

SR-1 (Tweezers-style)	10 Removers/Box	3
SR-3 (Scissors-style)	10 Removers/Box	3



21

## For more information

For more information about 3M Health Care Products:

- Contact your local 3M sales representative.
- Call the 3M Health Care Customer Helpline: 1 800 228 3957.
- Visit our 3M Health Care products web site: <http://www.3M.com/healthcare>



**Health Care  
3M Canada**

Post Office Box 5757  
London, Ontario N6A 4T1

1 800 563 2921 Ontario  
1 800 361 2650 Quebec  
1 800 361 4488 Atlantic Provinces  
1 800 661 8126 Western Canada

Printed in Canada  
© 3M 2000  
0007 MS12376 E

# 3M Skin Stapler Systems



*Innovation in skin closure*

**3M** *Innovation*