



## NUMERICAL INDEX

005.....	24	3130-27.....	8	9313-27.....	18
006.....	24	3130-30.....	8	9313-30.....	18
009.....	24	3200-27.....	9	9704-27.....	22
016.....	40	3200-30.....	9	9704-30.....	22
030.....	25	3300-27.....	9	9705-27.....	22
040.....	25	3300-30.....	9	9705-30.....	22
077.....	24	3600-24.....	40	9777-27.....	20
088.....	25	3600-27.....	40	9777-30.....	20
20.....	25	3610-24.....	40	Ready Room	
60.....	25	3610-27.....	40	1013-27-RR.....	32
100-27.....	15	3620-24.....	38	3600-27-RR.....	35
100-30.....	15	3620-27.....	38	3620-27-RR.....	35
200-27.....	14	3713-10.....	40	7020-RR.....	37
200-30.....	14	3713-15.....	40	8890-RR.....	36
300-27.....	15	3713-16.....	40	8890-RRFF.....	34
300-30.....	15	3770-10.....	40	81360-RR.....	36
1010-27.....	13	3770-15.....	40		
1010-30.....	13	3770-16.....	40	COLOR SELECTORS:	
1013-27.....	10	6190.....	13	In an effort to stay current, all color	
1013-30.....	10	6196.....	13	selections for individual products can be	
1020-27.....	12	6198.....	6	found online at: <a href="http://www.clinton-ind.com">www.clinton-ind.com</a>	
1020-30.....	12	8834.....	27		
1030-27.....	13	8844.....	26		
1030-30.....	13	8850.....	19		
1300-27.....	15	8860.....	31		
1300-30.....	15	8870.....	30		
3010-27.....	6	8890.....	29		
3010-30.....	6	9013-27.....	17		
3013-27.....	6	9013-30.....	17		
3013-30.....	6	9044-27.....	18		
3020-27.....	7	9044-30.....	18		
3020-30.....	7	9074-27.....	17		
3030-27.....	4	9074-30.....	17		
3030-30.....	4	9105-30.....	17		
3100-27.....	9	9133-27.....	16		
3100-30.....	9	9133-30.....	16		

## WHY CLINTON?

Buyers of medical equipment are faced with evaluating clinical needs and justifying expenditures on nearly every purchase. Now, Clinton makes it a whole lot easier on everybody with products that maximize performance and aesthetics while minimizing costs.

**Performance.** With nearly forty years in medical manufacturing, Clinton has met the evolving challenges of the industry with products that fit the needs of providers, while at the same time maintaining a track record of durability that boasts the best warranty in the business.

**Aesthetics.** Nobody likes outdated anything, so Clinton has changed its products over the years to reflect the advances in medical practice, design and contemporary decor, making products with appeal in both form and function.

**Cost.** Bells and whistles are great, but they also drive up costs. On the other hand, why pay more if you don't have to? Clinton products fill needs instead of empty pockets, providing medical professionals with just the right products—at just the right prices.



### IMPORTANT NOTE TO DISTRIBUTORS DOING BUSINESS IN CALIFORNIA:

Clinton Industries products drop shipping or to our knowledge being reshipped into California will carry appropriate Safe Harbor Warnings in compliance with California Proposition 65.

For more information on California Proposition 65, go to: [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)



## ALPHA SERIES TREATMENT TABLES

Clinton Alpha Series Treatment Tables are sleek, smart looking and affordable. From their smooth powder-coated frame to their durable steel construction and knit-backed upholstery, the Alpha Series tables represent one of the best values in their class.

*Alpha Series, Treatment Tables feature:*

- Gray, powder-coated, all-welded, steel frame
- Adjustable height levelers
- Pre-threaded hardware for easy assembly
- Paper dispenser included
- Rounded corner top
- 2" of firm foam padding (5 cm)
- Adjustable backrest standard on all models
- Choice of upholstery colors
- Fast, easy assembly with basic tools



### **Straight Line, Treatment Table with Shelving**

- Rigid all steel frame
- Adjustable backrest
- Gray, Laminate shelving
- 450 lbs. load capacity under normal use (204 kg)

3030-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

**3030-30** 30" width model (76 cm)

**Flat Top** See model 3300 on page 9



## ALPHA SERIES TREATMENT TABLES

### Straight Line, Treatment Table

- Rigid all steel frame
- Adjustable backrest
- 450 lbs. load capacity under normal use (204 kg)

3010-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

3010-30 30" width model (76 cm)

Flat Top See model 3100 on page 9



3010

### Straight Line, Treatment Table with Shelf and 2 Drawers

- Rigid all steel frame
- Adjustable backrest
- Gray, laminate shelf & drawer unit
- Easy-clean, one-piece, plastic drawers
- 450 lbs. load capacity under normal use (204 kg)

3013-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

3013-30 30" width model (76.2 cm)

Flat Top See model 3130 on page 8



3013

### Bariatric, Treatment Table

- Plywood top and welded steel frame
- 6-leg design with bolt-through assembly
- Leg levelers
- Gray, laminate full length shelf with H-braces
- 600 lbs. load capacity under normal use (272 kg)

6198	length	width	height
	72"	30"	31"
	183 cm	76 cm	79 cm



6198



### Straight Line, Treatment Table with Shelf

- Rigid, all steel frame
- Adjustable backrest
- Gray laminate shelf
- 450 lbs. load capacity under normal use (204 kg)

3020-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

3020-30

Flat Top See model 3200 on page 9

## ALPHA SERIES FLAT TOP TREATMENT TABLES

Clinton, Alpha Series, flat top, Treatment Tables combine classic steel construction and traditional upholstery with Clinton affordability to make them a very popular choice for an economical, durable treatment table.

Alpha Series, Flat Top, Treatment Tables feature:

- Gray, powder-coated, all-welded, steel frame
- Adjustable height levelers
- Pre-threaded hardware for easy assembly
- Paper dispenser included
- 2" of firm foam padding (5 cm)
- Choice of 20 upholstery colors
- Rounded corner top
- 5-year, limited warranty
- Fast, easy assembly with basic tools
- 500 lbs. load capacity under normal use (227 kg)

### Straight Line, Flat Top, Treatment Table with Shelf & Two Drawers

- Rigid, all steel frame
- Gray laminate shelf & drawer unit
- 450 lbs. load capacity under normal use (204 kg)

3130-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

3130-30	30" width model (76 cm)
---------	-------------------------



### Straight Line, Flat Top, Treatment Table

- Rigid all steel frame
- 450 lbs. load capacity under normal use (204 kg)

3100-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

3100-30	30" width model (76.19 cm)
---------	----------------------------



3100

### Straight Line, Flat Top, Treatment Table with Full Shelf

- Rigid all steel frame
- Gray, laminate shelf
- Easy-clean, one-piece, plastic drawers
- 450 lbs. load capacity under normal use (204 kg)

3200-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

3200-30	30" width model (76.19 cm)
---------	----------------------------



3200

### Straight Line, Flat Top, Treatment Table with Shelving

- Rigid all steel frame
- Gray, Laminate shelving
- 450 lbs. load capacity under normal use (204 kg)

3300-27	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.73 cm

3300-30	30" width model (76.19 cm)
---------	----------------------------



3300



3130

## CLASSIC SERIES TREATMENT TABLES

Clinton Classic Series Treatment Tables are some of our most popular tables. Tough, sturdy 2 1/4" hardwood legs support bolted-frame construction for the ultimate in strength and durability in a wood, treatment table.

*Classic Series, Treatment Tables feature:*

- Solid hardwood legs
- H-brace construction
- Legs have pre-inserted carriage bolts for easy assembly
- Rounded corner top
- 2" of firm foam padding (5 cm)
- Adjustable backrest standard on all models
- Available in Natural & Dark Cherry finish
- Choice of upholstery colors
- Fast, easy assembly with basic tools



### **Straight Line, Treatment Table with Shelf and Drawers**

- Hardwood legs & H-brace construction
- Adjustable backrest
- Laminate shelf & drawer unit
- 400 lbs. load capacity under normal use (181 kg)

1013-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

**1013-30** 30" width model (76 cm)

**Flat Top** See model 3100 on page 15



1013



1013

1013 table shown with optional paper dispenser



1020

**Straight Line, Treatment Table with Shelf**

- Hardwood legs & H-brace construction
- Adjustable backrest
- Full laminate shelf
- 400 lbs. load capacity under normal use (181 kg)

1020-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

1020-30 30" width model (76 cm)

Flat Top See model 200 on page 14



1020

**Straight Line, Treatment Table**

- Hardwood legs & H-brace construction
- Adjustable backrest
- 400 lbs. load capacity under normal use (181 kg)

1010-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

1010-30 30" (76 cm) width model

Flat Top See model 100 on page 15



1010

**Straight Line, Treatment Table with Shelving**

- Hardwood legs & H-brace construction
- Adjustable backrest
- Two-level laminate shelving
- 400 lbs. load capacity under normal use (181 kg)

1030-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

1030-30 30" (76 cm) width model

Flat Top See model 300 on page 15



1030

**Bariatric, Treatment Table with Shelf**

- Plywood top and solid hardwood legs
- 6-leg design with triple bolted corner legs
- Built-in middle leg levelers
- Steel reinforced frame
- Full length, laminate shelf with H-braces
- Pillow included
- 600 lbs. load capacity under normal use (272 kg)

6190	length	width	height
	78"	30"	31"
	198 cm	76 cm	79 cm

**Low Height, Bariatric, Treatment Table**

- Same as 6190 with 26" height
- 600 lbs. load capacity under normal use (272 kg)

6196	length	width	height
	78"	30"	26"
	198 cm	76 cm	66 cm



6190

## CLASSIC SERIES FLAT TOP TREATMENT TABLES

Clinton Classic Series Flat top Treatment Tables feature traditional wood construction and tough, durable upholstery with Clinton affordability to make them one of our most popular choices in a flat top treatment table.

Classic Series, Flat Top, Treatment Tables feature:

- Solid hardwood legs
- H-brace construction
- Legs have pre-inserted carriage bolts for fast easy assembly
- Seamless, rounded corner top
- 2" of firm foam padding (5 cm)
- Available in Natural and Dark Cherry finish
- Choice of 20 upholstery colors
- 5-year, limited warranty
- Fast, easy assembly with basic tools
- 400 lbs. load capacity under normal use (158.75 kg)



### Straight Line, Flat Top, Treatment Table with Shelf

- Hardwood legs & H-brace construction
- Full laminate shelf
- 400 lbs. load capacity under normal use (181 kg)

200-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

200-30	30" width model (76 cm)		
--------	-------------------------	--	--

### Straight Line, Flat Top, Treatment Table

- Hardwood legs & H-brace construction
- 400 lbs. load capacity under normal use (181 kg)

100-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

100-30	30" (76 cm) width model		
--------	-------------------------	--	--



100

### Straight Line, Flat Top, Treatment Table with Shelving

- Hardwood legs & H-brace construction
- Two-level laminate shelving
- 400 lbs. load capacity under normal use (181 kg)

300-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

300-30	30" width model (76 cm)		
--------	-------------------------	--	--



300

### Straight Line, Flat Top, Treatment Table with Shelf and Drawers

- Hardwood legs & H-brace construction
- Laminate shelf & drawer unit
- 400 lbs. load capacity under normal use (181 kg)

1300-24	length	width	height
	72"	24"	31"
	183 cm	61 cm	79 cm

1300-27	27" width model (69 cm)		
---------	-------------------------	--	--



1300





## STYLE LINE TREATMENT TABLES

Clinton Style Line, Treatment Tables feature easy-clean, flat panel, laminate surfaces and sleek, modern styling, making them a popular alternative to traditional table design.

Style Line, Treatment Tables feature:

- 1 1/8" panel legs with nylon glides (3 cm)
- Vertical spanner for extra support
- Bolt-through, pre-threaded hardware
- Easy-clean laminate surfaces
- Chip-resistant, PVC edges
- Rounded corner top
- 2" of firm foam padding (5 cm)
- Adjustable backrest standard on all models
- Available in your choice of laminate colors
- Choice of upholstery colors
- Fast, easy assembly with basic tools



### Treatment Table with Stool

- 3 drawers, storage shelf, adjustable backrest and matching stool that stores neatly under the shelf
- Stool is 17.75" wide x 17.75" deep (sits on floor) (45 cm x 45 cm)
- Self-closing drawers with metal sides and 75 lbs. capacity (34 kg)
- 400 lbs. load capacity under normal use (181 kg)

9133-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

9133-30 30" width model (76 cm)

Flat Top Option Use -088 after model number

### Cabinet Style, Treatment Table with Drawers & Doors

- Features 3 drawers, adjustable backrest and 1 storage compartment with 1 adjustable shelf
- Flat panel drawers with all steel sides
- Adjustable doors with concealed soft-close hinges
- 425 lbs. load capacity under normal use (193 kg)

9105-30	length	width	height
	72"	30"	31"
	183 cm	76 cm	79 cm

Flat Top Option Use -088 after model number



9105-30

### Treatment Table with Drawers & Shelf

- Two easy clean, plastic drawers, shelf and adjustable backrest
- 400 lbs. load capacity under normal use (181 kg)

9013-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

9013-30 30" width model (76 cm)

Flat Top Option Use -088 after model number



9013

### Cabinet Style, Treatment Table with 4 Doors

- Features 4 doors, adjustable backrest and 2 storage compartments each with an adjustable shelf
- Adjustable doors with concealed soft-close hinges
- 425 lbs. load capacity under normal use

9074-27	length	width	height
	72"	274"	31"
	183 cm	691 cm	79 cm

9074-30 30" width model (76 cm)

Flat Top Option Use -088 after model number



9074



9133

## STYLE LINE TREATMENT TABLES

### Cabinet Style, Treatment Table with 4 Doors

- Features an adjustable backrest, 2 storage compartments, 2 sliding doors and an adjustable shelf in each compartment
- Adjustable doors with concealed soft-close hinges
- 425 lbs. load capacity under normal use



9044-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

9044-30 30" width model (76 cm)

Flat Top Option Use -088 after model number



9313

### Flat Top Style Line Straight Line Treatment Table with Shelf and Two Drawers

- Steel reinforced, pre-assembled plywood top frame
- Coordinating laminate two drawer unit and shelf
- 1 1/8" panel legs with nylon glides (2.85 cm)
- Vertical spanner for extra support
- Bolt-through, pre-threaded hardware
- Fast, easy assembly with basic tools
- 400 lbs. load capacity under normal use (181.43 kg)

9313-27	length	width	height
	72"	27"	31"
	183 cm	69 cm	79 cm

30" width model (76 cm)

Flat Top Option Use -088 after model number



8850

### Panel Leg, Space-Saver, Treatment Table

- Panel leg design
- Easy-clean, all laminate construction
- Pull-out, laminate, leg rest extends length of the table
- Adjustable pneumatic backrest
- 400 lbs. load capacity under normal use (181 kg)

8850	length	width	height
	58"/72" extended	28 1/2"	31"
	147 cm/183 cm extended	72 cm	81 cm