

---

---

## MATERIAL SAFETY DATA SHEET

**ProAdvantage Latex Textured Powder Free Exam Gloves**  
P359101, P359102, P359103, P359104, P359105

### Physical Data – Composition

<b>Natural Solid Rubber</b>	<b>75.8%</b>
<b>Sulfur</b>	<b>1.0%</b>
<b>Zinc Accelerator</b>	<b>1.2%</b>
<b>Others</b>	<b>22.0%</b>

### Fire and Explosions

Stability: Stable  
Extinguishing Media: Foam, Carbon Dioxide, Water  
Flash Point: N/A

### Normal Precautions

Ventilation: Not Necessary Under Normal Conditions  
Spillage Handling: N/A  
Eyes: N/A  
Storage: Store at Room Temperature

### Precaution

Do not reuse gloves. Good quality exam gloves provide an excellent biological barrier. They are not intended for use as a chemical barrier.

### Storage

Store in original packing in a cool, dry and well ventilated area, away from dust, sunlight, moisture, X-ray and excessive heat about 100° F (37° C)

---

---