



Industrial/Emergency

First Aid Kit

The Pro Advantage® first aid kit is perfect for a wide variety of workplaces. Whether serving construction sites, medical facilities, manufacturing facilities or offices, this kit has been designed to contain the critical supplies necessary for first aid emergencies.

OSHA regulations require the availability of adequate first aid supplies for all employees at the workplace. Made of durable plastic, this kit features quality products meeting both OSHA regulations and ANSI guidelines.



FEATURES

- Meets federal OSHA regulation 1910.151b (states may vary)
- Meets ANSI Z308.1-2003 when minimum fill is maintained
- Serves workplace up to 25 people
- Trilingual labels in English, Spanish and French
- Contains first aid essentials to prepare for workplace accidents
- Back contents label for retail environment
- Includes First Aid Guide
- Wall mountable
- Made in the USA

CONTENTS – 158 PIECES

- 1 ea Plastic case with dividers, 9" x 8-3/8" x 2-1/2"
- 100 ea 1" x 3" Adhesive plastic bandages*
- 3 ea Knuckle fabric bandages
- 3 ea Fingertip fabric bandages
- 1 set 40" x 40" x 56" Triangular sling/bandage with 2 safety pins*
- 4 ea 3" x 3" Gauze pad, 2 2-packs*
- 5 ea 2" x 3" Non-stick pads
- 1 ea 5" x 9" Trauma pad*
- 1 ea 3" Conforming gauze roll bandage
- 1 ea 2" Conforming gauze roll bandage
- 2 ea Sterile eye pads
- 15 ea Antiseptic cleansing wipes (sting free)*
- 3 ea Ammonia inhalant ampoules
- 6 ea First aid/burn cream packs, 0.9 gm*
- 1 ea 4" x 5" Instant cold compress
- 1 ea 1/2" x 10 yd. First aid tape roll*
- 1 ea 4-1/2" Scissors, nickel-plated
- 1 ea 4" Plastic tweezers, one time use
- 4 ea Exam quality vinyl gloves, 2 pairs*
- 1 ea Eye wash, 1/2 oz.
- 1 ea First Aid Guide

*Required Minimum Fill

Item #	Description
P440025	First Aid Kit, 25 Person, Plastic Case