

Plain and Frosted

# Microscope Slides and Cover Glass



## MICROSCOPE SLIDE FEATURES

- Guaranteed pre-cleaned
- Thin, flat and resistant to corrosion or fogging
- Made from the highest quality white Swiss glass
- Individually selected for uniform quality
- Precision ground edges on all four sides
- Frosted slides feature a sandblasted double frosted end for ease of use
- Slides measure 25mm x 75mm x 1mm

## COVER GLASS FEATURES

- Free from ragged edges, oily film, striae, and finger prints
- Measures 22mm x 22mm

Item #	Description	Quantity
P460025	Plain Microscope Slides	72/box
P460125	Frosted Microscope Slides	72/box
P460122	Cover Glass (0.13 to 0.17mm thick)	1 ounce/box
P460222	Cover Glass (0.17 to 0.25mm thick)	1 ounce/box