

Aneroid

Sphygmomanometers



Deluxe Aneroid Sphygmomanometer

Offering the most features, the Pro Advantage Deluxe Aneroid Sphygmomanometer includes a precision crafted manometer, cuff size marking system with extended graduation marks, filter protected bulb, dip molded inflation bladder for heavy use, chrome-plated brass air release valve and a zippered carry case. Lifetime calibration warranty, 2 year system warranty. Latex-free.

P548200 Aneroid, Adult, Black P548210 Aneroid, Small Adult, Black P548215 Aneroid, Large Adult, Black Premium Pocket Aneroid Sphygmomanometer

The Pro Advantage Premium Pocket Aneroid Sphygmomanometer is a professional inflation system featuring a manometer, a cuff size marking system, chrome-plated brass air release valve, dip molded inflation bladder and bulb for heavy use and a zippered carry case. Lifetime calibration warranty, 2 year system warranty. Latex-free.

P548340 Aneroid, Adult, Black P548330 Aneroid, Small Adult, Black P548360 Aneroid, Large Adult, Black P548380 Aneroid, Thigh, Black



Standard Pocket Aneroid Sphygmomanometer

The Pro Advantage Standard Pocket Aneroid Sphygmomanometer is an economical inflation system featuring a manometer, chrome-plated deflation valve, nylon cuff, inflation bladder and bulb and a zippered carry case. 10 year calibration warranty, 1 year system warranty. Latex-free

P548440 Aneroid, Adult, Navy Blue P548445 Aneroid, Small Adult, Navy Blue P548450 Aneroid, Large Adult, Navy Blue

Aneroid

Sphygmomanometers







Cuff & Bladder



Inflation System

Sphygmomanometer Accessories

Accessories for Pro Advantage Sphygmomanometers. Available choices are a pocket gauge, cuffs and bladders, and inflation systems. Gauge features a black enamel finish. Cuff and bladder features a nylon cuff with marking system and a dip molded bladder. Inflation system features a nylon cuff with marking system, inflation bladder, filter protected bulb and valve. Two year warranty; Pocket Gauge has a lifetime calibration warranty. Latex-free.

P548900	Pocket Gauge, Black	P549515	Inflation System, Adult, Black
P549220	Cuff & Bladder, Adult, 2-Tube, Black	P549510	Inflation System, Small Adult, Black
P549240	Cuff & Bladder, Large Adult, 2-Tube, Black	P549525	Inflation System, Large Adult, Black
P549260	Cuff & Bladder, Thigh, 2-Tube, Black	P549530	Inflation System, Thigh, Black

Distributed By:

