

DD-701

Rechargeable Display Doppler

Newman Medical

Newman Medical is proud to offer a full line of Dopplers with the performance and durability you expect to meet your application - obstetrical or vascular. The **DigiDop 701 rechargeable, display** Doppler has unique features not found in other Doppler systems such as Digital Optimized Processing (DOP™). Each DigiDop is designed with the sensitivity you need while also designed to be tough and last through years of everyday office use. All our products are **Made in the USA in Colorado**. Each Newman Medical Doppler and probe must pass a series of comprehensive tests to assure quality and the durability of each product.



Features:

- ✦ **Best Sound** - Digitally Optimized Processing – DOP™ produces up to 6x signal of other Dopplers
- ✦ **Best Warranty** – 5 year \$0 cost warranty – **ALL** Parts. **ALL** Labor.
- ✦ **Best Support & Service** – Knowledgeable Staff & Guaranteed Service
- ✦ Tough and durable polycarbonate and stainless steel case
- ✦ 'Unbreakable' magnetic probe holder
- ✦ All probes are interchangeable
- ✦ Each main unit comes with choice of one probe
 - ✦ 2MHz, 2MHz WP, 3MHz, 3MHz WP, 5MHz, 8MHz
- ✦ User replaceable cord
- ✦ Powered by 3 AA batteries offer extended battery life (900 mins)
- ✦ Dimensions (h x w x l): 150 x 65 x 35 mm (6.0 x 2.5 x 1.25 in)
- ✦ Weight: 340 grams (12 ounces)



Newman Medical
5350 Vivian Street | Unit C | Arvada | CO | 80002
Toll-free 800-267-5549 | Fax-888-202-3627
www.newman-medical.com

NM0004B