

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

EPD PACKAGE 2

EPIDURAL POSITIONING DEVICE & STAND

STANDARD FEATURES:

- Load capacity: 600 lbs.
- 20.5"W x 30"D x 50-60"H (Assembled)
- White powder coat base
- Black upholstery
- Includes 1 Box (50) Sani-Cover® Fitted Disposable Face Rest Covers
- Epidural Positioning Device (EPD)
- 180° Adjustable face rest
- 6 Position adjustable arm rest
- Arm rest size: 16.5"W x 10"D
- Arm rest adjustment: 32-43" from the floor
- Interchangeable torso cushions (1 flat and 1 elevated)
- Torso cushion size: 8"W x 8"H
- Torso cushion adjustment: 5.5-10"
- EPD Stand (LDS)
- Positive locking clamp for EPD
- 7 Position adjustable foot rest
- Foot rest size: 20"W x 12"D
- Foot rest adjustment: 10.5-19.5" from the floor
- Medical grade profile locking casters

EPD Package 2
(Epidural Positioning
Device + Stand)
with flat torso cushion



Elevated torso
cushion



FRC 30-50



Patents:
US 6,758,447,B2
and D455,495

STANDARD PACKAGE:

☐ EPD Package 2 Epidural Positioning Device & EPD Stand
QTY:



EPD - Epidural
Positioning Device



LDS - EPD Stand