



## WALL & FREESTANDING CABINETS

Cabinets are a great way to securely store your AED in public.

Cabinets increase device visibility and allow for easy access in times of emergency. All cabinets are International Building Code certified and comfortably fit all models of Physio-Control AEDs. Alarm and strobe light features help ensure your AEDs are protected.

### Color Choices



White



Stainless Steel



OSHA Orange



### Surface Mount (Projects 7" from Wall)

Cat. Number	Color	Material	Size	Edge	Alarm	Strobe	Fire Resistant
11998-000310	White	Carbon Steel	17.38" x 17.38" x 7"	Square	•	•	
11220-000076	Stainless Steel	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•		
11220-000079	White	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•		
11220-000082	OSHA Orange	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•		
11220-000083	White	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•	•	
11220-000084	Stainless Steel	Carbon Steel	14" x 14" x 6.75"	Rolled	•	•	
11220-000086	OSHA Orange	Aluminum	17.38" x 17.38" x 7"	Square	•	•	
11220-000087	OSHA Orange	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•	•	



### Semi-Recessed (Projects 3" from Wall)

---

Cat. Number	Color	Material	Size	Edge	Alarm	Strobe	Fire Resistant
11998-000292	White	Carbon Steel	14" x 14" x 6.75"	Rolled	•		
11210-000026	White	Carbon Steel	14" x 14" x 6.75"	Rolled	•		•
11220-000077	Stainless Steel	Carbon Steel	14" x 14" x 6.75"	Rolled	•		



### Recessed Mount (Projects 1.5" from Wall)

---

Cat. Number	Color	Material	Size	Edge	Alarm	Strobe	Fire Resistant
11210-000027	White	Carbon Steel	14" x 14" x 6.75"	Square	•		•
11220-000078	Stainless Steel	Carbon Steel	14" x 14" x 6.75"	Square	•		
11220-000088	OSHA Orange	Aluminum	14" x 14" x 6.75"	Square	•	•	
11220-000089	OSHA Orange	Carbon Steel	14" x 14" x 6.75"	Square	•		
11220-000090	OSHA Orange	Carbon Steel	14" x 14" x 6.75"	Square	•	•	
11998-000293	White	Carbon Steel	14" x 14" x 6.75"	Square	•		



### Free Standing

---

Cat. Number	Color	Material	Size	Edge	Alarm	Strobe	Fire Resistant
11210-000028	White	Carbon Steel	17.38" x 17.38" x 7"	Square	•		
11210-000029	Gray	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•		

## AED LOCATION WALL SIGNS

AED Location Wall Signs alert responders to the location of an AED.

AED Location Wall Signs come in various designs and models to match the environment while meeting local AED recommendations. Lighter weight materials are less obtrusive and easier to hang while brighter designs improve visibility. Help decrease response times by signaling the location of your AED.



Tent AED Location Sign



Flat AED Location Sign



T-mount AED Location Sign

Cat. Number	Style	Material	Size Wide x Tall	Thickness	Weight	Depth from Wall	Width across Wall
11998-000332	Tent w/ Logo	White Styrene	11" x 8"	0.079"	4 oz	4" deep	11"
11998-000330	Flat w/ Logo	White Styrene	8.5" x 11"	0.02"	1 oz	Flush	8.5"
11998-000331	T-mount w/ Logo	White PVC	8" x 10"	0.125"	4 oz	8" deep	1" wide lip
11998-000333	T-mount w/o Logo	White PV	8" x 10"	0.125"	4 oz	8" deep	1" wide lip

---

DISTRIBUTED BY:



Contact your local Red Cross representative for more information.

