

# Wide Variety of Test Strips Provides Medical Professionals Flexibility

The CardioChek Plus Lipid+eGLU® Smart Bundle Pack™ Includes:  
PTS Panels® lipid panel test strips and eGLU test strips



## Three-Analyte



### Lipid Panel Test Strip

Measures: Total Cholesterol, HDL Cholesterol & Triglycerides  
Calculates: LDL Cholesterol, Total Cholesterol/HDL Ratio, Non-HDL Cholesterol



### CHOL+HDL+GLU Panel Test Strip

Measures: Total Cholesterol, HDL Cholesterol & Glucose  
Calculates: Total Cholesterol/HDL Cholesterol Ratio, Non-HDL Cholesterol



### METABOLIC CHEMISTRY Panel Test Strip

Measures: HDL Cholesterol, Triglycerides & Glucose

## Dual-Analyte



### CHOL+GLU Panel Test Strip

Measures: Total Cholesterol & Glucose



### CHOL+HDL Panel Test Strip

Measures: Total Cholesterol & HDL Cholesterol  
Calculates: Total Cholesterol/HDL Cholesterol Ratio

## Single-Analyte



### TOTAL CHOLESTEROL Test Strip

Measures: Total Cholesterol



### HDL CHOLESTEROL Test Strip

Measures: HDL Cholesterol



### TRIGLYCERIDES Test Strip

Measures: Triglycerides



### GLUCOSE Test Strip

Measures: Glucose



### eGLU Test Strip

Measures: Electrochemical Glucose

Single-analyte test strips available as special order include Creatinine and Ketone.



**PTS Diagnostics**  
7736 Zionsville Road, Indianapolis, IN 46268 USA  
Office: +1 (317) 870-5610 FAX: +1 (317) 870-5608  
Toll-Free in the US: 1 (877) 870-5610