



Clinton takes power tables to new heights with a complete line of Power, Hi-Lo Tables designed and engineered to fit the needs of physicians, technicians and therapists.

Clinically proven in the field, Clinton Power Tables, are among the best in their class. With quality construction, dependable smooth, operation and enough power to satisfy nearly every need, Clinton Power Tables are practical, and best of all, affordable.

Welded steel construction makes them tough and two inches of firm foam padding topped with premium, knit-backed vinyl makes patient comfort a priority in all our models. Choose your Power, Hi-Lo Table by selecting a load capacity to fit your needs in either Shrouded or Open Base models. Clinton also offers Open Base Power Imaging Tables or Power Ultrasound and Imaging Trendelenburg Tables to meet specific needs.

Make Clinton your one source for affordable, convenient Power, Hi-Lo Tables — and plug into the power of Clinton.

INDEX

Shrouded, Power Tables2-7
 Open Base, Power Tables.....8-15
 Power XL Tables16-17
 Open Base, Power, Imaging Table18-19
 E-Series Power Tables.....20-21
 Power Glide, Treatment Tables 22-23
 Open-Base, Power Changing Tables 24-25
 Power, Ultrasound, Trendelenburg Tables26-29
 Power, Imaging, Trendelenburg Tables30-33
 Power Table Options 34-35
 Ready Room™36-37

IMPORTANT NOTE TO DISTRIBUTORS DOING BUSINESS IN CALIFORNIA:

Clinton Industries products drop shipping or to our knowledge being reshipped into California will carry appropriate Safe Harbor Warnings in compliance with California Proposition 65.

For more information on California Proposition 65, go to: www.P65Warnings.ca.gov



Clinton power tables are among the most practical and affordable power exam tables on the market. They offer smooth dependable power, stylish design and dependable service.

Clockwise from left to right: 81360–Shrouded Power Table with Adjustable Backrest, Drop Section & Stirrups; 81330–Shrouded Power Table with Adjustable backrest & Drop Section; 80079–Multi-Use, Power Imaging, Trendelenburg Table with Adjustable Backrest, Drop Section & Stirrups; 80330 Open Base Power Imaging Table with Adjustable Backrest & Drop Section.



POWER 500

Shrouded, Power Tables feature:

- Ultra-strong, uni-frame constructed base and top frame provides a rigid patient platform
- All welded steel base and top frame
- 3/4" plywood top for extra support (1.9 cm)
- Textured ABS plastic encloses lift mechanism and makes cleaning easier
- Smooth electric height adjustment moves straight up and down
- Pneumatic, backrest cylinder is internally mounted on the frame for a cleaner appearance
- Hardened bolts and self lubricating pivot points for years of service
- Dependable 110 volt motor
- Hands-free foot control
- Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights may help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism
- Many different top configurations and options available
- Ships fully assembled with no upcharge
- 500 lbs. combined load and force capacity under normal use (227 kg)

POWER 500, Shrouded, Power Table with One Piece Top

- Optional side rails, locking casters & more (pages 30 & 31)

81100	length	width	height range
overall	72"	27"	21"-34"
	183 cm	69 cm	53 cm - 86 cm

POWER 500, Shrouded, Power Table with Adjustable Backrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Optional side rails, locking casters & more (pages 30 & 31)

81200	length	width	height range
overall	72"	27"	21"-34"
	183 cm	69 cm	53 cm - 86 cm
backrest	38"	27"	
	97 cm	69 cm	



81200



81100

⚠ WARNING: Products featured on this page can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P6SWarnings.ca.gov.



Shown with optional, solid drop section

POWER 500, Shrouded, Power Table with Adjustable Backrest & Drop Section

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (specify when ordering)
- Optional side rails, locking casters & more (pages 30 & 31)



81330

81330	length	width	height range
overall	76"	27"	21"-34"
	193 cm	69 cm	53 cm - 86 cm
backrest	38"	27"	
	97 cm	69 cm	
drop section	13"	27"	
	33 cm	69 cm	

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P6SWarnings.ca.gov.

SHROUDED, POWER TABLES

POWER 500, Shrouded, Space Saver Power Table with Drop Section

- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- Optional side rails, locking casters & more (*pages 30 & 31*)

81220	length	width	height range
overall	76"	27"	21"-34"
	193 cm	69 cm	53 cm – 86 cm
drop section	13"	27"	
	33 cm	69 cm	

POWER 500, Shrouded, Power Table with Adjustable, Casting Table with ClintonClean™ Leg Rest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- ClintonClean™ Solid Plastic Surface leg rest is adjustable from either side and raises up to 30° from flat position
- Optional side rails, locking casters & more (*pages 30 & 31*)

81350	length	width	height range
overall	76"	27"	21"-34"
	193 cm	69 cm	53 cm – 86 cm
backrest	38"	27"	
	97 cm	69 cm	
leg rest	13"	26"	
	33 cm	66 cm	

POWER 500, Shrouded, Power Table with Stirrups, Adjustable Backrest & Footrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and raises up to 30° from flat position
- Retractable, 2-position, adjustable stirrups
- Optional side rails, locking casters & more (*pages 30 & 31*)

81360	length	width	height range
overall	76"	27"	21"-34"
	193 cm	69 cm	53 cm – 86 cm
backrest	38"	27"	
	97 cm	69 cm	
footrest	13"	19"	
	33 cm	48 cm	



⚠ WARNING: Products featured on this page can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

SHROUDED, EXTRA-WIDE, BARIATRIC, POWER TABLES

POWER 600 BARIATRIC

POWER 600 Bariatric, Shrouded, Extra-Wide Straight Power Table

- All welded, uni-frame steel base and top frame
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*pages 30 & 31*)

84108-34	length	width	height range
overall	72"	34"	21"-34"
	183 cm	86 cm	53 cm – 86 cm

84108-40 40" (101.60 cm) width model

POWER 600 Bariatric, Shrouded, Extra-Wide, Power Table with Adjustable Backrest

- Dual cylinder, pneumatic backrest, adjustable from either side and raises up to 80° from flat
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*pages 30 & 31*)

84208-34	length	width	height range
overall	72"	34"	21"-34"
	183 cm	86 cm	53 cm – 86 cm
backrest	38"	34"	
	97 cm	86 cm	

84208-40 40" (101.60 cm) width model

POWER 600 Bariatric, Shrouded, Extra-Wide, Power Table with Adjustable Backrest & Drop Section

- Dual cylinder pneumatic backrest, adjustable from either side and raises up to 80° from flat
- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*pages 30 & 31*)

84438-34	length	width	height range
overall	76"	34"	21"-34"
	193 cm	86 cm	53 cm – 86 cm
backrest	38"	34"	
	97 cm	86 cm	
drop section	13"	34"	
	33 cm	86 cm	

84438-40 40" (101.60 cm) width model



⚠ WARNING: Products featured on this page can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

POWER 500

POWER 500, Open Base, Power Tables feature:

- Ultra-strong uni-frame constructed base and top frame provides a rigid patient platform
- 3/4" plywood top for extra support (1.9 cm)
- All welded steel base, top frame and scissors lift
- Smooth electric height adjustment moves straight up and down
- Pneumatic, backrest cylinder are internally mounted on the frame for a cleaner appearance
- Hardened bolts and self lubricating pivot points for years of service
- Dependable 110 volt motor
- Hands-free foot control
- Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights may help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism
- Many different top configurations and options available
- Ships fully assembled with no upcharge
- 500 lbs. combined load and force capacity under normal use (227 kg)

POWER 500, Open Base, Power Table with Stirrups, Adjustable Backrest & Footrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Retractable, 2-position, adjustable stirrups
- Optional side rails, locking casters & more (pages 30 & 31)

80360	length	width	height range
overall	76" 193 cm	27" 69 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	27" 69 cm	
footrest	13" 33 cm	19" 48 cm	



POWER 500, Open Base, Power Table with One Piece Top

- Optional side rails, locking casters & more (pages 30 & 31)

80100	length	width	height range
	72" 183 cm	27" 69 cm	18"-35" 46 cm - 89 cm



POWER 500, Open Base, Power Table with Adjustable Backrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Optional side rails, locking casters & more (pages 30 & 31)

80200	length	width	height range
overall	76" 193 cm	27" 69 cm	18"-35" 46 cm - 89 cm
Backrest	38" 97 cm	27" 69 cm	



POWER 500, Open Base, Power Table with Adjustable Backrest & Drop Section

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section, adjustable from either side and lowers up to 80° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (specify when ordering)
- Optional side rails, locking casters & more (pages 30 & 31)

80330	length	width	height range
overall	76" 193 cm	27" 69 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	27" 69 cm	
drop section	13" 33 cm	27" 69 cm	



⚠ WARNING: Products featured on this page can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

⚠ WARNING: Products featured on this page can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

OPEN BASE, POWER TABLES

POWER 500, Space Saver, Power Table with Drop Section

- Drop section is adjustable from either side and lowers up to 80° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (*pages 30 & 31*)

80220	length	width	height range
overall	76"	27"	18"-35"
	193 cm	69 cm	46 cm - 89 cm
drop section	13"	27"	
	33 cm	69 cm	



POWER 500, Power Table, Adjustable Casting Table with ClintonClean™ Surface Leg Rest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position
- ClintonClean™ Solid Plastic Surface leg rest is adjustable from either side and raises up to 30° from flat position
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (*pages 30 & 31*)

80350	length	width	height range
overall	76"	27"	18"-35"
	193 cm	69 cm	46 cm - 89 cm
backrest	38"	27"	
	97 cm	69 cm	
leg rest	13"	26"	
	33 cm	66 cm	



⚠ WARNING: Products featured on this page can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

OPEN BASE, EXTRA-WIDE, BARIATRIC, POWER TABLES

POWER 600 BARIATRIC

Open Base, Extra-Wide, Bariatric, Straight Power Table

- All welded, uni-frame steel base and top frame
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*pages 30 & 31*)

84100-34	length	width	height range
overall	72"	34"	18"-35"
	183 cm	86 cm	46 cm - 89 cm

84100-40 Same as above with 40" width (102 cm)



Open Base, Power Table with Adjustable Backrest

- Dual cylinder, pneumatic backrest, adjustable from either side and raises up to 80° from flat
- All welded, uni-frame steel base and top frame
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*pages 30 & 31*)

84200-34	length	width	height range
overall	72"	34"	18"-35"
	183 cm	86 cm	46 cm - 89 cm
backrest	38"	34"	
	97 cm	86 cm	

84200-40 Same as above with 40" widths (102 cm)



Open Base, Power Table with Adjustable Backrest & Drop Section

- Dual cylinder pneumatic backrest, adjustable from either side and raises up to 80° from flat
- Drop section is adjustable from either side and lowers up to 80° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*pages 30 & 31*)

84430-34	length	width	height range
overall	76"	34"	18"-35"
	193 cm	86 cm	46 cm - 89 cm
backrest	38"	34"	
	97 cm	86 cm	
drop section	13"	34"	
	33 cm	86 cm	

84430-40 Same as above with 40" widths (102 cm)



OPEN BASE, POWER XL TABLES

CLINTON OPEN BASE & SHROUDED
POWERXL
 Extra Width, Extra Capacity, Extra Affordable

Power XL Tables feature:

- 600 lbs. (272 kg) load capacity under normal use
- 5-year limited warranty on frame and electric motor drive
- Smooth, electric, raise and lower height adjustment
- Pneumatic backrest cylinder is mounted internally
- Heavy duty, steel frame is fabricated in the USA by Clinton Industries
- Dependable scissors-lift features an H-leg design for extra stability
- 30" Rounded corner 3/4" plywood top with 2" of ultra firm padding covered by heavy weight upholstery
- Harden steel bolts and self lubricating pivot bearings for years of dependable service
- Low height may help achieve ADA tax credit
- Optional side rails, locking casters & more (pages 30 & 31)



80199



80299



80399

Open Base, Power XL Table with One Piece Top

- Top is 76" long x 30" wide

80199	length	width	height range
overall	72" 183 cm	30" 76 cm	18"-35" 46 cm - 89 cm

Open Base, Power XL Table with Adjustable Backrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position

80299	length	width	height range
overall	76" 183 cm	30" 76 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	30" 76 cm	

Open Base, Power XL Table with Adjust. Backrest and Drop Section

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position

85309	length	width	height range
overall	72" 183 cm	30" 76 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	30" 76 cm	
footrest	13" 33 cm	30" 76 cm	

SHROUDED, POWER XL TABLES

Shrouded, Power XL Table with One Piece Top

- Top is 76" long x 30" wide
- Textured ABS plastic encloses lift mechanism and makes cleaning easier

80199	length	width	height range
overall	72" 183 cm	30" 76 cm	18"-35" 46 cm - 89 cm



81199

Shrouded, Power XL Table with Adjustable Backrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Textured ABS plastic encloses lift mechanism and makes cleaning easier

80299	length	width	height range
overall	76" 183 cm	30" 76 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	30" 76 cm	



81299

Shrouded, Power XL Table with Adjust. Backrest and Drop Section

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Textured ABS plastic encloses lift mechanism and makes cleaning easier

85309	length	width	height range
overall	72" 183 cm	30" 76 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	30" 76 cm	
footrest	13" 33 cm	30" 76 cm	



81399

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

OPEN BASE, POWER IMAGING TABLES

Open Base, Power Imaging Tables feature:

- Ultra-strong uni-frame constructed base and top frame provides a rigid patient platform
- 3/4" plywood top for extra support (1.9 cm)
- All welded steel base, top frame and scissors lift
- Smooth electric height adjustment moves straight up and down
- Pneumatic, backrest cylinder are internally mounted on the frame for a cleaner appearance
- Hardened bolts and self lubricating pivot points for years of service
- Built-in hinged window drop. 10" x 10"
- Metal slide bar supports window drop
- Dependable 110 volt motor
- Hands-free foot control
- Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights may help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism that does not need to roll on the floor
- Many different top configurations and options available
- 5-year limited warranty on frame and electrical components
- Ships fully assembled with no upcharge

Open Base, Multi-Use, Power Imaging Table with Stirrups

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Retractable, 2-position, adjustable stirrups
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (pages 30 & 31)

85309	length	width	height range
overall	76" 193 cm	27" 69 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	27" 69 cm	
footrest	13" 33 cm	19" 48 cm	



85309

Flat position with window drop engaged

Open Base, Power Imaging Table with Window Drop

- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (pages 30 & 31)

85100	length	width	height range
overall	72" 183 cm	27" 69 cm	18"-35" 46 cm - 89 cm



85100

Open Base, Power Imaging Table with Adjustable Backrest and Window Drop

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- All welded steel base, top frame and scissors lift
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (pages 30 & 31)

85200	length	width	height range
overall	72" 183 cm	27" 69 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	27" 69 cm	



85200

Power Imaging Table with 3 Section Top

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (pages 30 & 31)

85300	length	width	height range
overall	76" 193 cm	27" 69 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	27" 69 cm	
footrest	13" 33 cm	27" 69 cm	



85300

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

⚠ WARNING: Products featured on this page can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.



E-Series Power Tables feature:

- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 1.5" (4 cm) ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options (page 27)
- Made in the USA by Clinton Industries
- Backed by a 5-year limited warranty
- See options on page 31
- 500 lbs. (227 kg) load capacity under normal use

Open Base, Power Table with One Piece Top

86100	length	width	height range
	72"	27"	18"-35"
	183 cm	69 cm	46 cm - 89 cm



86100

Open Base, Power Table with Adjustable Backrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position

86200	length	width	height range
overall	72"	27"	18"-35"
	183 cm	69 cm	46 cm - 89 cm
Backrest	38"	27"	
	97 cm	69 cm	



86200

Open Base, Power Table with Adjustable Backrest & Drop Section

- Drop section, adjustable from either side and lowers up to 80° from flat position

86220	length	width	height range
overall	76"	27"	18"-35"
	193 cm	69 cm	46 cm - 89 cm
drop section	13"	27"	
	33 cm	69 cm	



86220

Open Base, Power Table with Adjustable Backrest & Footrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position

86300	length	width	height range
overall	76"	27"	18"-35"
	193 cm	69 cm	46 cm - 89 cm
backrest	38"	27"	
	97 cm	69 cm	
footrest	13"	19"	
	33 cm	48 cm	



86300

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

power-glide™
Hi-Lo Treatment Tables by Clinton

Power-Glide Tables feature:

- Smooth-glide, electric power height adjustment
- Ultra heavy-duty base with solid, one-piece, steel ends
- Easy-clean, ABS plastic, center shroud hides lift components and provides a clean look
- Rigid, one-piece, laser-cut steel top frame
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Fast easy motor access via ABS side panels
- Hands-free low voltage foot control
- Seamless, rounded corner, upholstered top
- 2" ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries
- 500 lbs. (227 kg) load capacity under normal use
- 3-year limited factory warranty from Clinton Industries



Power-Glide, Hi-Lo, Treatment Table with One Piece Top

- Features an easy-clean, one-piece top

87100	length	width	height range
overall	72" 183 cm	27" 69 cm	18.5"-34" 47 cm - 86 cm

Power-Glide, Hi-Lo, Treatment Table with Adjustable Backrest

- Features easy-clean, two-piece top with adjustable backrest
- Pneumatic backrest is adjustable from either side and raises to a maximum of approximately 80°

87200	length	width	height range
overall	72" 183 cm	27" 69 cm	18.5"-34" 47 cm - 86 cm
Backrest	38" 97 cm	27" 69 cm	



Power-Glide, Hi-Lo, Treatment Table with Flat Top and Drop Section

- Features a flat top and friction-lock, drop section
- Drop section is adjustable from either side and raises to approximately 75° degrees from flat.

87220	length	width	height range
overall	76" 193 cm	27" 69 cm	18.5"-34" 47 cm - 86 cm
drop section	13" 33 cm	27" 69 cm	



Power-Glide, Hi-Lo, Treatment Table with Adjustable Backrest & Drop Section

- Features an adjustable backrest and friction-lock, drop section
- Pneumatic backrest is adjustable from either side and raises to a maximum of approximately 80°
- Drop section is adjustable from either side and raises to approximately 75° degrees from flat.

87330	length	width	height range
overall	76" 193 cm	27" 69 cm	18.5"-34" 47 cm - 86 cm
Backrest	38" 97 cm	27" 69 cm	
drop section	13" 33 cm	27" 69 cm	



Power-Glide, Hi-Lo, Treatment Table with Adjustable Backrest, Footrest and Stirrups

- Features an adjustable backrest, friction-lock, footrest and stirrups
- Pneumatic backrest is adjustable from either side and raises to a maximum of approximately 80°
- Footrest is adjustable from either side and raises to approximately 75° from flat
- Retractable, stirrups swing out to a maximum of 35°

87360	length	width	height range
overall	76" 193 cm	27" 69 cm	18.5"-34" 47 cm - 86 cm
Backrest	38" 97 cm	27" 69 cm	
drop section	13" 33 cm	18" 46 cm	

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.



Only Clinton has the manufacturing know-how and the focus to engineer and bring to market a changing table that meets needs and saves money. Building on our reputation for dependable and affordable medical equipment, Clinton's Open Base, Power, Changing Table is easy to operate and comfortable for patients.



Open Base, Power Changing Tables feature:

- Easy-clean, one-piece top
- Two lock-in-place, padded side rails
- Four, smooth-roll, locking, 3" casters
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 2" (5 cm) ultra-firm padding
- Four, locking casters
- Low height may help achieve ADA tax credit
- Two Lengths: 72" and 60"
- Made in the USA by Clinton Industries
- 400 lbs. (181 kg) load capacity under normal use

⚠️ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.



80110-60 one side down



80110-60 both sides up



80110-72 one side down




80110-72 both sides down

60" Open Base, Power, Changing Table

80110-60	length	width	height range
overall	60"	27"	18" – 35"
	152 cm	69 cm	46 cm – 89 cm

087-6 EXTRA-HEIGHT OPTION

- Height extension to accommodate Hoyer Lift
- Adds 6.25" to the height of the table from top pad to floor



72" Open Base, Power, Changing Table

80110-72	length	width	height range
overall	72"	27"	18" – 35"
	183 cm	69 cm	46 cm – 89 cm

⚠️ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

Power, Ultrasound, Trendelenburg Tables feature:

- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.1
- Industry leading 5-year, limited warranty
- Standard models: 500 lbs load capacity under normal use (227 kg)
- Extra-wide models: 600 lbs load capacity under normal use (272 kg)

General, Flat Top, Ultrasound Table

- One piece flat top with firm foam padding
- Optional side rails, locking casters & more (pages 30 & 31)

80061	length	width	height range
	72"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
	183 cm	69 cm	55 cm – 95 cm

80061-X Extra-wide: 34" (86 cm)



80061



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

! WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

General Ultrasound Table with Three-Section Top

- Three-piece top with pneumatic adjustable backrest and manual footrest
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Optional side rails, locking casters & more (pages 30 & 31)

80063	length	width	height range
overall	76"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
	193 cm	68.58 cm	55 cm – 95 cm
backrest	38"	27"	
	97 cm	68.58 cm	
footrest	13"	27"	
	33 cm	69 cm	

80063-X Extra-wide: 34" (86 cm)



80063

POWER
Backrest & Footrest Available

! WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

Multi-Use, Ultrasound Table with Stirrups

- Three-piece top with pneumatic adjustable backrest, manual footrest, and stirrups
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Optional side rails, locking casters & more (pages 30 & 31)

80069	length	width	height range
overall	76"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
	193 cm	68.58 cm	55 cm – 95 cm
backrest	38"	27"	
	97 cm	69 cm	
footrest	12"	18.5"	
	30 cm	47 cm	

80069-X Extra-wide: 34" (86 cm)



80069

POWER
Backrest & Footrest Available

! WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

POWER, ULTRASOUND,
TRENDLENBURG TABLES

**General, Ultrasound Table with
Adjustable Backrest**

- Two-piece top with easy-to-adjust, pneumatic backrest
- Optional side rails, locking casters & more
(pages 30 & 31)

80062

overall	72"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
	183 cm	69 cm	55 cm – 95 cm

backrest	38"	27"	POWER Backrest Available
	97 cm	69 cm	

80062-X Extra-wide: 34" (86 cm)



⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

Power Imaging, Trendelenburg Tables feature:

- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.1
- Industry leading 5-year, limited warranty
- Standard models: 500 lbs load capacity under normal use (227 kg)
- Extra-wide models: 600 lbs load capacity under normal use (272 kg)

Flat Top, Power Imaging Table with Drop Window

- One piece flat top with drop window
- Optional side rails, locking casters & more (pages 30 & 31)

80071	length	width	height range
overall	72" 183 cm	27" 69 cm	21 1/2" – 37 1/2" 55 cm – 95 cm
window	10" 25 cm	10" 25 cm	

80071-X Extra-wide: 34" (86 cm)



80061



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

Power Imaging Table with Fowler Back and Drop Window

- Two-piece top with pneumatic adjustable backrest and window drop
- Optional side rails, locking casters & more (pages 30 & 31)

80072	length	width	height range
overall	72" 183 cm	27" 69 cm	21 1/2" – 37 1/2" 55 cm – 95 cm
backrest	38" 97 cm	27" 69 cm	
window	10" 25 cm	10" 25 cm	

80072-X Extra-wide: 34" (86 cm)



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

Power Imaging Table with Three-Section Top and Drop Window

- Three-piece top with pneumatic adjustable backrest, manual footrest and window drop
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Optional side rails, locking casters & more (pages 30 & 31)

80073	length	width	height range
overall	76" 193 cm	27" 68.58 cm	21 1/2" – 37 1/2" 54.61 cm – 95.25 cm
backrest	38" 96.52 cm	27" 68.58 cm	
footrest	12" 30.48 cm	27" 68.58 cm	
window	10" 25.4 cm	10" 25.4 cm	

80073-X Extra-wide: 34" (86 cm)



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters



80072

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.



80073

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

POWER, IMAGING,
TRENDELENBURG TABLES

Multi-Use, Power Imaging Table with Stirrups and Drop Window

- Three-section table with pneumatic adjustable backrest, manual footrest, stirrups and window drop
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Optional side rails, locking casters & more (pages 30 & 31)

80079	length	width	height range
overall	76"	27"	21 1/2" – 37 1/2"
	193 cm	69 cm	55 cm – 95 cm
backrest	38"	27"	
	97 cm	69 cm	
footrest	12"	18 1/2"	
	30 cm	47 cm	
window	10"	10"	
	25 cm	25 cm	



80079-X Extra-wide: 34" (86 cm)



⚠️ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

POWER TABLE OPTIONS

Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop hook and loop fasteners
- 2 straps per set

080	length	width
	92"	2"
	233.67 cm	5 cm

Side Rails (pair)

- 3 height position
- Plunger height adjustment

082 Standard / 098 Trendelenburg

Power Backrest

(Trendelenburg Tables only)

- Smooth power control raises backrest to 80° and lowers to flat

011

Power Footrest

(Trendelenburg Tables only)

- Smooth power control lowers footrest to 30° and raises to flat

090

Hand-Held Height Control

(Not available on E-Series, or Power-Glide Tables)

- Alternative to foot control
- Stores easily under top when not in use

084

Locking Casters

(Not available on E-Series Tables)

- Heavy duty, rubber wheel, easy-roll casters

087

30" Width Top

(Trendelenburg Tables only)

037



Paper Dispenser

(Standard Tables & Power-Glide Option)

- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030

Paper Cutter

(Standard Tables & Power-Glide Option)

- Clear plastic with elastic straps and buckle

040

Paper Dispenser/Cutter Combo

((Standard Tables & Power-Glide Option)

- Paper dispenser and paper cutter ordered as one option

034



ONLY OPTIONS AVAILABLE FOR E-SERIES POWER TABLES

Nose Hole

- Nose hole in drop section of 86220 or 86300

009



Paper Dispenser

- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030



Paper Cutter

- Clear plastic with elastic straps and buckle

040

Paper Dispenser/Cutter Combo

- Paper dispenser and paper cutter ordered as one option

034



E-Series
POWER TABLES

WARNING: See specific table product for California Proposition 65 Warning.

WARNING: See specific table product for California Proposition 65 Warning.



One Complete Exam Room, One Order Number, One Price

A complete office, well equipped and ready for patients. Clinton's Power Ready Room, or E-Series Ready Room packages offer big advantages in function, one source convenience and at a pricepoint that fits your budget.

81360-RR

(no model substitutions, your choice of Fashion Finish cabinetry and upholstery colors)

- 81360 in your choice of upholstery color
- 8348 Wall Cabinet in your choice of Fashion Finish
- 8648 Base Cabinets with sink in your choice of Fashion Finish
- 48P your choice of Postform Countertop
- 2135 in your choice of upholstery color
- C-40B Premium, Padded, Side Chair in your choice of upholstery color
- GW-3000 Single, Stainless Steel* Glove Box holder
- TR-20S 20 Quart, Medium, Round, Stainless Steel*, Waste Receptacle

* not MRI compatible

⚠️ WARNING: One or more of the products in this grouping can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.



86300-RR

(no model substitutions, your choice of Fashion Finish cabinetry and upholstery colors)

- 86300 in your choice of upholstery color
- 8348 Wall Cabinet in your choice of Fashion Finishes
- 8648 Base Cabinets with sink in your choice of Fashion Finishes
- 48P your choice of Postform Countertop
- 2135 in your choice of upholstery color
- C-40B Premium, Padded, Side Chair in your choice of upholstery color
- GS-3000 Single, Stainless Steel* Glove Box holder
- TR-20S 20 Quart, Medium, Round, Stainless Steel* Waste Receptacle


* not MRI compatible

⚠️ WARNING: One or more of the products in this grouping can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

For more information go to www.P65Warnings.ca.gov.

fashionfinish

COUNTERTOP COLORS





800-441-9131 • 717-848-3519

Toll-Free Fax 866-522-1233

Fax 717-843-5871

www.clinton-ind.com