




Triage[®]

Rapid OXY | BUP | MDMA PANEL


Rapid PCP Single TEST

Get additional testing flexibility with Triage Rapid visual-read test cassettes for Oxycodone, Buprenorphine, MDMA, and PCP.


- CLIA waived
- Dip and read or pipette method.
- Quick results in approximately 3 to 5 minutes.
- Results are stable for up to 6 hours.
- Room temperature (15°C to 30°C).



Analyte	Test Abbreviation	Level (ng/mL)
Buprenorphine	BUP	12.5
MDMA ((+/-) 3,4-methylenedioxy-methamphetamine)	MDMA	1000
Oxycodone	OXY	100
Phencyclidine (phencyclidine CI)	PCP	25



Dip Method



Pipette Method

Cat. #	Description	Kit Size
16-20348	Triage Rapid OXY BUP MDMA Panel	25 Tests
16-20349	Triage Rapid PCP Single Test	25 Tests
10-ZNC-000	Detectabuse [®] Universal Negative Control	1 x 5.0 mL
10-ZPCOX-220	Detectabuse Universal Positive Control	1 x 5.0 mL

