

QUICKVUE[®]

One-Step hCG
Urine and Combo TEST

Two Simple Tests You Can Depend On For The Early Detection of Pregnancy

Rapid Diagnosis

- Results in 5 minutes (serum), 3 minutes (urine).
- Positive may become visible sooner.

Excellent Performance

- >99% sensitivity; specificity; accuracy.
- 25 mIU/mL sensitivity in serum or urine.

Convenient

- Room temperature storage.
- 18-month shelf life.

Easy to Use


- Simple, one-step procedure uses only 3 drops of sample.
- Two-color endpoint for easy-to-read and interpret results.

Reliable

- Highly sensitive assay with built-in controls.

Lot-to-Lot Consistency

- Highly controlled and monitored manufacturing environment.



QUIDEL[®]

For more information on QuickVue One-Step hCG Urine and Combo test kits, visit our website at www.quidel.com, or call your Quidel Account Manager at 1.800.874.1517 (U.S. only) or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • quidel.com

Combo Procedure:


Urine Procedure:



To order: Call your Quidel Account Manager at 1.800.874.1517.
Outside the U.S. call 858.552.1100.

QuickVue hCG Urine – 25 Test Kit:
Catalog # 20109
hCG Control Set (Urine):
Catalog # 00272

QuickVue hCG Combo – 50 Test Kit:
Catalog #20110
Serum hCG Control Set:
Catalog #00281



QUIDEL®

For more information on QuickVue One-Step hCG Urine and Combo test kits, visit our website at www.quidel.com, or call your Quidel Account Manager at 1.800.874.1517 (U.S. only) or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • quidel.com

Distributor Information: