



TMX-428

High-quality, low-maintenance treadmill

High-quality, low maintenance

MEDICAL TREADMILL

This heavy-duty medical treadmill is specifically designed to be interfaced with and controlled by the **SCHILLER** Stress Test system and a wide variety of cardiac and pulmonary stress testing systems.

The **TMX-428** is distinguished by its smooth, quiet, and reliable performance, while its design prioritizes the absolute safety and comfort of the patient. The treadmill is self-calibrating, virtually maintenance-free, and features a low-profile deck and a 0.1 MPH minimum speed for enhanced safety.


FEATURES

- Smooth, quiet, zero-start, 4HP (6HP peak) hi-torque drive
- 500 lb. patient capacity
- Smooth, quiet linear actuator elevation system to 25%
- Low-to-the ground, cushioned deck for patient comfort and safety
- 22" x 63" (56cm x 160cm) running belt and safety handrails
- Belt alignment system eliminates off-center running belts
- Exceptionally accurate, self-calibrating speed and elevation
- Prominent emergency stop button and patient safety tether
- Minimum maintenance



The **TMX-428** Treadmill is specifically engineered to meet the rigorous demands of medical motion devices. Its quiet operation makes it an ideal choice for both hospital and office settings, ensuring a smooth and unobtrusive experience for patients and healthcare providers.

Technical Specifications



Base dimension	84" x 36" (L x W)
Weight	480 lbs.
Max. Patient weight	500 lbs. / 227 kg
Safe, Low-profile deck	7" (18 cm) from floor
Running surface	22" x 63" (56 x 160 cm)
Self-aligning running belt	Yes
Cushioned and reversible running deck	Yes
All steel chassis	Yes
Powder Coat finish	Yes
Emergency stop and patient tether	Yes
Space Requirements	33" x 78.5" (84 x 199 cm)
Stress Test Interface	RS-232

Drive System

Heavy-duty, 4 HP (6 HP peak) drive

Power supply | 208-240V, 50/60 Hz, dedicated 15A service
110V, 50/60 Hz, dedicated 20A service

Speed Range

Self-calibrating | Yes

Elevation Range | 0-25%, Self-calibrating

Other Options

Emergency stop button left or right side

Adjustable front handrail

SCHILLER

A M E R I C A S



Call Us:
+1(786) 485-0620



Visit us at:
www.schillerus.com



10903 NW 33rd Street
Doral FL, 33172, USA.

