seca 403 and seca 402

Carts for mobile support of seca baby scales

new



- Custom-made indentations on the surface hold seca baby scales firmly in place.*
- Ideal working height of 31.5 inches / 80 cm.
- Integrated drawer and shelf offer extra storage space so you always have everything within reach (seca 403).
- Optional available paper roll holder helps you keep the work area clean (seca 402).



seca 403 and seca 402:

The perfect supplement to seca baby scales.

They were specially developed for safe and hygienic routine use in hospitals and medical practices. Customized indentations on the surface accommodate the scale's base and ensure that the scale does not slip out of place. Two practical brakes keep the cart stable so babies can be safely weighed. And you can push or pull the cart by its handle to move the baby scale quickly to the place where it's needed.

For work that's simple, ergonomic and hygienic.

Mobile and safe use of baby scales is of prime importance in pediatrics. So now the seca baby scales* can be securely positioned on specially developed carts. The ergonomic design of the carts, with a work surface over 31 inches, ensures safe work processes. Hygienic standards are easier to maintain, thanks to the seamless, easy-to-clean surface.

The seca 403 offers practical storage options.

The expanded variant of the seca cart has practical storage space that keeps documents and materials close at hand. Extra space is provided by a shelf and an integrated drawer. Safe measuring and weighing anywhere – with the seca 402.

In the basic model the seca 402 cart offers all the advantages of safe, mobile and hygienic work with seca baby scales. A practical and simple to install paper roll holder is optionally available.

The optional paper roll holder on the seca 402 makes work even easier and more hygienic.





the seca baby scale anchored on the cart.



Four bumpers on the wheels protect the cart and scale from collisions with walls and doors.

seca 403 and seca 402

Technical Data

seca 403:

- Dimensions (WxHxD): 28.2 x 31.5 x 17.1" / 715 x 800 x 435 mm
- Weight: 60.6 lbs / 27.5 kg
- Max. load: 220 lbs / 100 kg
- 4 bumpers
- 2 brakes
- Integrated drawer and shelf

seca 402:

- Dimensions (WxHxD): 28.2 x 31.5 x 17.1" / 715 x 800 x 435 mm
- Weight: 25.4 lbs / 11.5 kg
- Max. load: 220 lbs / 100 kg
- 4 bumpers
- 2 brakes
- Optional: paper roll holder
- * Both carts offer an individual surface for following seca baby scales: seca 374, seca 334, seca 354, seca 728 and seca 727.

North America West • seca corp. • 13601 Benson Avenue • Chino, CA 91710 USA North America East • seca corp. • 7240 Parkway Drive, Suite 120 • Hanover, MD 21076 USA Toll free +1 800 542 7322 • Fax +1 888 705 7397 • info.us@seca.com • www.seca.com