

# EMR validated column scale with eye-level display

**seca**

Precision for health

## seca 787

- + Medical-grade, robust scale for years of use with high patient volume
- + Digital measuring rod with an extra long measuring range
- + Two-line display shows height and weight
- + EMR validated: Send results directly into the EMR
- + Damping system and Auto-HOLD guarantee quick and precise weighing
- + Display swappable to left or right side

 EMR  
Integration



seca 787

# EMR validated column scale with eye-level display

## Advanced EMR validated weighing technology

The seca 787 fully meets the demands of professional EMR (electronic medical record) integrated measuring and weighing. It boasts a tilt-proof full metal platform and a stable single piece column that is at eye level and won't give way when operating the keypad. It is a quality every patient will notice, as soon they set foot on the platform. The scale is easy to move by staff due to the large rubber coated castors and a wide wheelbase.



## Absolute sturdy construction

With its absolute sturdiness, the seca 787 gives patients (especially those affected by obesity) a greater sense of security when stepping on the platform. The digital measuring rod is firmly attached to the column and provides an extra long measuring range. A generous two-line display shows weight and height simultaneously and automatically calculates patients' BMI. The damping system and auto-HOLD function enable quick and precise measurements.

## seca integration solutions:

How workflows can clearly be improved.

The seca 787 can be integrated into your information system via the network. At the push of a button, measured values can be sent directly to the EMR, where they are immediately available to all users. Extended with an optional USB barcode scanner, both patient and user are identified and the measurement results are always correctly assigned. This saves time, increases process reliability and eliminates transmission errors.

Further information on seca integration solutions at: [integration.seca.com](http://integration.seca.com)



Capacity	550 lbs / 250 kg
Increments	0.2 lbs / 100 g
Dimensions (WxHxD)	16.1 x 53.4 x 22.7 inch / 41.0 x 135.6 x 57.6 cm
Platform (WxHxD)	13.8 x 16.5 x 2.3 inch / 35.0 x 41.8 x 5.9 cm
Measuring range	4 – 90 inch / 10 – 230 cm (digital measuring range 48 – 90 inch / 122 – 230 cm)
Weight	28.0 lbs / 12.7 kg

Power supply	Batteries, power adapter (optional)
Data transmission	LAN or Wi-Fi via interface module seca 452 (optional)
Functions	TARE, HOLD, Auto-HOLD, Damping (motion filter), lbs/kg switch-over lock, automatic BMI, auto-off
Electromagnetic radiation resistance	According to IEC 60601-1-2

We reserve the right to make modifications.

seca north america west  
seca corp.  
13601 Benson Avenue  
Chino – CA 91710 – USA

P +1 800 542 7322  
F +1 888 705 7397  
E [info.us@seca.com](mailto:info.us@seca.com)

Precision for health  
[seca.com](http://seca.com)

201178716019