

seca robusta 813

High capacity electronic flat scales.



- Solid steel frame, up to 440 lbs capacity.
- Non-slip rubber base for safe standing during measuring.
- Large LCD figures for easy reading.
- Step-off function for weighing without switching scales.
- From the global market leader for medical weighing and measuring.

seca robusta 813:

High capacity electronic flat scales for personal use.

The seca robusta 813 has an extra wide weighing platform that easily hold up to a 440 lbs capacity. The tough rubber coating takes the heaviest challenges lightly and the flat construction makes stepping on very easy.



A scale of sturdiness - with a touch of elegance.

Thanks to its solid steel frame, the seca robusta 813 can withstand 440 lbs and features large LCD figures for easy reading. Furthermore, its extra-wide surface and flat shape facilitate trouble-free weighing of heavyweight patients.

Pleasant and comfortable weighing.

All platforms are not created equal. For this weighing platform seca has chosen a high-quality, two-component nubby material which is soft like rubber and pleasant to the touch. But don't be fooled - this material is also tough. It fulfils two requirements at the same time: pleasant weighing and an easy-to-clean surface that stands up to all commercially available cleaning and disinfecting products.

Paying attention to details.

Weight and girth are no problem on the seca robusta 813. You don't even have to bend over to switch it on, thanks to the unique step-off function. An additional problem with many scales is the display. Thanks to the 1.1 inch-high LCD figures, it is extremely easy to read the measurements. You see, a scale from seca makes life easier by paying attention to details.



High-quality, two-component nubby material feels as soft as rubber.



Big LCD figures. For easy reading.

seca robusta 813

Technical data

- Capacity: 440 lbs / 200 kg
- Division: 0.2 lbs / 100 g
- Dimensions (WxHxD):
17 x 1.9 x 14.7" / 433 x 47 x 373 mm
- Weight: 6.3 lbs / 2.9 kg
- Functions: auto-HOLD, step-off, lbs/kg switch-over, tip-on automatic switch-on, automatic switch-off
- Parts to be delivered:
4 AAA batteries,
4 extension feet for carpets