




The first scale that knows your name





seca Scale-up Line

Developed for electronic medical records

We're not only concerned with the perfect integration of measured data into your electronic medical record, we also deliver the ideal measuring station. That way, you get a completely coordinated integration solution from one source for safer and more efficient workflows.

World innovation: the ID-Display™ is directly connected to the electronic medical record

The ID-Display™ lets you see which patient was recorded with the scanner to prevent any confusion. That means more control for you when transferring to the Electronic Medical Records (EMR) and significantly improved safety for users and patients.

Absolute confidence

- + The display shows names or IDs of users and patients
- + Display-controlled data transfer
- + Voice response for patient guidance aids personnel
- + Rapid measurement process within seconds thanks to the innovative ultrasound technology
- + Elegant design featuring extreme robustness thanks to a steel handrail directly anchored to the cast iron base
- + Large platform made of seca Bearclaw Glass® virtually break-resistant glass for extreme durability and hygiene




EMR Integration We've turned intricate into easy

Save time, eliminate transfer errors and achieve a zero-incident environment. Transfer height and weight directly into the Electronic Medical Records (EMR) and revolutionize your workflows. seca scales and measuring stations can be integrated into all EMR systems. Feel free to contact us.



Large platform made of seca Bearclaw Glass® virtually break-resistant for absolute durability and hygiene.



Convenient transport with extremely robust rubber wheels.



The optional panda motif makes taking children's measurements lots of fun.



Footprint LED lights show the patient the correct position for height measurement.



Always the highest quality equipment for your needs



EMR-validated flat scale with ID-Display™ and stable glass platform

+ Flat scale seca 650

+ Capacity: 360 kg / 800 lbs
+ Graduation: 50 g / 0.1 lb

Order code: **SECA 650**



EMR-validated handrail scale with ID-Display™

+ Flat scale seca 654
+ Handrail seca 455

+ Capacity: 360 kg / 800 lbs
+ Graduation: 50 g / 0.1 lb

Order code: **ON SHM INT NN**



EMR-validated ultrasonic measuring station with ID-Display™

+ Flat scale seca 654
+ Ultrasonic height measurement seca 257

+ Capacity: 360 kg / 800 lbs
+ Graduation: 50 g / 0.1 lb
+ Measuring range: 60 – 220 cm / 24 – 87"
+ Graduation of measuring range: 1 mm / 1/8"

Order code: **ON SMM IUT NN**



EMR-validated ultrasonic measuring station with ID-Display™ and handrail

+ Flat scale seca 654
+ Handrail seca 455
+ Ultrasonic height measurement seca 257

+ Capacity: 360 kg / 800 lbs
+ Graduation: 50 g / 0.1 lb
+ Measuring range: 100 – 220 cm / 40 – 87"
+ Graduation of measuring range: 1 mm / 1/8"

Order code: **ON SHM IUT NN**

seca headquarters

+49 40 20 00 00 0

seca gmbh & co. kg
Hammer Steindamm 3–25
Germany

seca middle east

+971 564 410 292

Dubai
United Arab Emirates

seca asia pacific

+60 173 711 973

Kuala Lumpur
Malaysia

seca.com/scaleup