



A PROGRESSIVE LEGACY OF

PRECISION PERFORMANCE



Ophthalmic Product Offering

40 Year Heritage Of Premium Performance

Quality Innovative Needle Designs

A Sharp Alternative

SURGICAL
SPECIALTIES
CORPORATION™



CORNEAL TRANSPLANT

Penetrating Keratoplasty (PK)

NYLON

ULTRAGLIDE

A9000N	10-0	30cm	DOUBLE ARMED	5.51mm,	1/2 circle	160 curvature
A9003N	10-0	30cm	DOUBLE ARMED	6.15mm,	3/8 circle	135 curvature
A2675N	10-0	30cm	DOUBLE ARMED	4.88mm,	bi-curve	90/50 curvature
A9060N	10-0	10cm	SINGLE ARMED	6.15mm,	3/8 circle	135 curvature
A9061N	10-0	10cm	SINGLE ARMED	5.51mm,	1/2 circle	160 curvature
A9001N	10-0	15cm	DOUBLE ARMED	5.51mm,	1/2 circle	160 curvature

NYLON

SIDE CUTTING LANCET

A7756N	10-0	30cm	DOUBLE ARMED	5.51mm,	1/2 circle	150 curvature
A7770N	10-0	30cm	DOUBLE ARMED	6.15mm,	3/8 circle	135 curvature
A7718N	10-0	30cm	DOUBLE ARMED	6.60mm,	3/8 circle	131 curvature
A2533N	10-0	30cm	SINGLE ARMED	6.15mm,	3/8 circle	135 curvature
A2532N	10-0	30cm	SINGLE ARMED	6.60 mm,	3/8 circle	131 curvature
A2534N	10-0	15cm	SINGLE ARMED	5.51mm,	1/2 circle	150 curvature




BLEPHAROPLASTY

POLYSN

SIDE CUTTING LANCET

G506N	5-0	45cm	SINGLE ARMED	16mm,	1/2 circle	175 curvature
G493N	5-0	45cm	SINGLE ARMED	13mm,	3/8 circle	135 curvature
G492N	6-0	45cm	SINGLE ARMED	13mm,	3/8 circle	157 curvature
G570N	6-0	45cm	SINGLE ARMED	8.51mm,	1/4 circle	109 curvature



IOL SUTURE FIXATION

POLYPROPYLENE

SIDE CUTTING LANCET

J2558N	9-0	15cm	DOUBLE ARMED	14.99mm,	1/8 circle	50 curvature
J2559N	10-0	7.5cm	DOUBLE ARMED	14.99mm,	1/8 circle	50 curvature



GLAUCOMA FILTERING SURGEY

Trabeculectomy and Shunt

POLYSN

SIDE CUTTING LANCET

GA-546N	7-0	45cm	DOUBLE ARMED	6.6mm,	3/8 circle	131 curvature
GA-566N	7-0	30cm	DOUBLE ARMED	6.60mm,	3/8 circle	131 curvature
GA-548N	8-0	30cm	DOUBLE ARMED	6.60mm,	3/8 circle	131 curvature

POLYSN

PRECISION REVERSE CUTTING

GA-8760N	9-0	13cm	SINGLE ARMED	4.40mm,	3/8 circle	145 curvature
----------	-----	------	--------------	---------	------------	---------------

POLYSN

ULTRAGLIDE

GA-448N	10-0	30cm	DOUBLE ARMED	5.50mm,	1/2 circle	160 curvature
---------	------	------	--------------	---------	------------	---------------




VITRECTOMY

NYLON

ULTRAGLIDE

A9000N	10-0	30cm	DOUBLE ARMED	5.51mm,	1/2 circle	160 curvature
A9003N	10-0	30cm	DOUBLE ARMED	6.15mm,	3/8 circle	135 curvature
A2675N	10-0	30cm	DOUBLE ARMED	4.88mm,	bi-curve	90/50 curvature
A9060N	10-0	10cm	SINGLE ARMED	6.15mm,	3/8 circle	135 curvature



GENERAL USAGE

SILK

SIDE CUTTING LANCET

D7733N	7-0	45cm	DOUBLE ARMED	6.60mm,	3/8 circle	131 curvature
DA2550N	8-0	45cm	DOUBLE ARMED	6.60mm,	3/8 circle	131 curvature
D7729N	8-0	45cm	SINGLE ARMED	14mm,	3/8 circle	131 curvature