










Nylon Sutures

Side Cutting Lancet Needles

Needle	Description	Material	Length	Size			
				10-0	9-0	8-0	
BSL5  TG6-C	Bicurve Length 4.88mm Curvature 90/50	Radius 1.52/2.79mm Chord 3.71mm Wire .15mm	Black Monofilament	6"/15cm	A2687N 7799G ●		
				12"/30cm	A2522N ●		
BSL5  TG6-C	Bicurve Length 4.88mm Curvature 90/50	Radius 1.52/2.79mm Chord 3.71mm Wire .15mm	Black Monofilament	8"/20cm	A2524N ●		
				12"/30cm	A2525N 7799G ●		
DSL6  TG140-6	3/8 Circle Length 6.15mm Curvature 135	Radius 2.54mm Chord 4.83mm Wire .15mm	Black Monofilament	6"/15cm	A2683N ●		
				12"/30cm	A2533N ●		
DSL6  TG140-6	3/8 Circle Length 6.15mm Curvature 135	Radius 2.54mm Chord 4.83mm Wire .15mm	Black Monofilament	12"/30cm	A7770N 7770G ●		
DSL7  TG140-8	3/8 Circle Length 6.60mm Curvature 131	Radius 2.79mm Chord 5.26mm Wire .20mm	Black Monofilament	12"/30cm	A2532N ●		
DSL7  TG140-8	3/8 Circle Length 6.60mm Curvature 131	Radius 2.79mm Chord 5.26mm Wire .20mm	Black Monofilament	12"/30cm	A7718N 7718G ●	A7717N 7717G ●	
HSL5  TG160-6	1/2 Circle Length 5.51mm Curvature 150	Radius 2.03mm Chord 4.07mm Wire .15mm	Black Monofilament	6"/15cm	A2534N ●	A2681N ●	
				12"/30cm	A2510N ●		
HSL5  TG160-6	1/2 Circle Length 5.51mm Curvature 150	Radius 2.03mm Chord 4.07mm Wire .15mm	Black Monofilament	6"/15cm	A7757N 7757G ●		
				12"/30cm	A7756N 7756G N2756K ●	A7760N 7760G ●	
HSL7  TG160-6	1/2 Circle Length 7.04mm Curvature 175	Radius 2.29mm Chord 4.71mm Wire .15mm	Black Monofilament	12"/30cm	A2519N ●		


Product numbers are color coded respectively as Sharpoint™ (dark blue), Ethicon™ (black), and Syneture™ (green).
B or N denotes a box of 12; D denotes a box of 36.

Nylon Sutures

UltraGlide™ Needles

Needle	Description	Material	Length	Size		
				11-0	10-0	9-0
BGL5  CSB-6	Bicurve Length 4.88mm Curvature 90/50 Radius 1.52/2.79mm Chord 3.71mm Wire .15mm	Black Monofilament	12"/30cm		A2675N 9007G ●	
DGL6  CS140-6 SE-140-6	3/8 Circle Length 6.15mm Curvature 135 Radius 2.54mm Chord 4.82mm Wire .15mm	Black Monofilament	4"/10cm		A9060N 9060G ●	
DGL6  CS140-6 SE-140-6	3/8 Circle Length 6.15mm Curvature 135 Radius 2.54mm Chord 4.82mm Wire .15mm	Black Monofilament	12"/30cm		A9003N 9003G ●	
HGL5  CS160-6 SE-160-6	1/2 Circle Length 5.51mm Curvature 160 Radius 2.03mm Chord 4.06mm Wire .15mm	Black Monofilament	4"/10cm		A9061N 9061G ●	
HGL5  CS160-6 SE-160-6	1/2 Circle Length 5.51mm Curvature 160 Radius 2.03mm Chord 4.06mm Wire .15mm	Black Monofilament	6"/15cm	A9016N 9016G ●	A9001N 9001G ●	
			12"/30cm		A9000N 9000G ●	
HGL7  CS175-8	1/2 Circle Length 6.99mm Curvature 175 Radius 2.29mm Chord 4.57mm Wire .15mm	Black Monofilament	12"/30cm		A9032N 9032G ●	

UltraGlide™ Needles

Needle	Description	Material	Length	Size		
				10-0	9-0	8-0
VVL4 	1/4 Circle Length 4mm 100 Microns Curvature 109	Black Monofilament	5"/13cm		A0819N ● ●	

Product numbers are color coded respectively as SharpPoint™ (dark blue), Ethicon™ (black), and Syneture™ (green).
B or N denotes a box of 12; D denotes a box of 36.