






PolySyn™ (Polyglycolic Acid) Sutures

Precision Reverse Cutting Needles (cont'd)

Needle	Material	Length	Size						
			8-0	7-0	6-0	5-0	4-0	3-0	
HSM9 9mm 1/2 Circle  P-2 P-21	Undyed Braided	18"/45cm					G503N J503G SL1644G	G504N J504G SL1645G	
HSM16 16mm 1/2 Circle  PS-4 P-24	Undyed Braided	18"/45cm					G506N J506G SL5632G	G507N J507G SL5633G	

Side Cutting Lancet Needles

Needle	Description	Material	Length	Size				
				8-0	7-0	6-0	5-0	
VSL7  TG100-8 SE-100-8	1/4 Circle Length 6.60mm Curvature 94	Radius 3.94mm Chord 5.89mm Wire .20mm	Violet Braided	18"/45cm			G544N J544G L1750K	
VCL8  S-14 SS-14	1/4 Circle Length 8.51mm Curvature 112	Radius 4.36mm Chord 7.24mm Wire .33mm	Undyed Braided	18"/45cm			G670N J670G L1670	G671N J671G L1669
	1/4 Circle Length 8.51mm Curvature 109	Radius 4.32mm Chord 7.28mm Wire .34mm	Violet Braided	8"/20cm			G590N J590G L1740K	
				18"/45cm			G570N J570G ●	
VSP8  S-29, S-24 SD-1, SS-24	1/4 Circle Length 8.51mm Curvature 112	Radius 4.36mm Chord 7.24mm Wire .33mm	Violet Braided	12"/30cm			G556N J556G L1446K	
	1/4 Circle Length 8.51mm Curvature 109	Radius 4.32mm Chord 7.28mm Wire .34mm		18"/45cm			G555N J555G ●	
				12"/30cm			G552N J552G L1742K	G553N J553G L1785K

Product numbers are color coded respectively as SharpPoint™ (dark blue), Ethicon™ (black), and Syneture™ (green).
 B or N denotes a box of 12; D denotes a box of 36.