



Sharp^{point}TM PLUS
Surgical Suture

Choice. Collaboration. Value.

Sharp^{point}TM PLUS Surgical Suture Catalog

SURGICAL
SPECI⁽ALTIES
CORPORATIONTM

For the last 40 years, Surgical Specialties Corporation has been trusted for quality, and we continue to strive to meet the needs and expectations of surgical professionals worldwide.

SharpPoint™ PLUS is the preferred global brand for physicians' wound and surgical closure needs where choice, collaboration and value matter.

The portfolio features a number of specialty sutures designed to meet the exacting needs of plastic, dermatologic, orthopedic and general surgeons.

SharpPoint™ PLUS Surgical Sutures are offered in a variety of absorbable and non-absorbable materials and attached to a wide variety of needle options.

The value and quality of SharpPoint™ PLUS sutures are augmented by a team of sales representatives dedicated to supporting the needs of our customers.

Table of Contents


How to Use Our Catalog	1
Materials Chart	2
A Guide to Surgical Sutures.....	3-4
Needle Anatomy and Terminology	5
Needle Shapes and General Uses	6
Needle Points and Cutting Edges	7
Needle Charts.....	8-10

Absorbable Sutures

Plain Gut	11-12
Chromic Gut.....	13-15
PolySyn FA™ (polyglycolic acid).....	16
Monoderm™ (glycolide and ε-caprolactone)	17-18
PolySyn™ (polyglycolic acid).....	19-24
PDO (polydioxanone).....	25

Non-Absorbable Sutures

Silk	26-29
Nylon™	30-33
Polypropylene	34-35
Polyviolene™ (polyester) - Uncoated	36



Choice. Collaboration. Value.

Introduction

How to Use Our Catalog

The pages in this catalog show the following information for each product code. Some products feature an expanded chart with detailed measurements (not shown in this sample) to ensure proper needle selection.

Needle Description - Contains an illustration of the needle. The comparable Ethicon™ needle code is black, and the comparable Syneture™ needle code is green. This section also lists the needle length (decimal 2 places included if not equal to zero) and shape.

Needle Geometry - Type of needle with tip geometry

Material - Color and style of suture

Reverse Cutting

Length - Length in inches/centimeters

Size Grid - Indicates the U.S.P. size for each product code listed.

Needle	Material	Length	Size		
			6-0	5-0	4-0
DSM12 12mm 3/8 Circle C-2 C-1	Black Braided Silk	18"/45cm	D1639N 1639G S1172	S1173	736G S1174

•• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10"

— Sharpoint™ PLUS product code
 — Ethicon™ comparable product code
 — Syneture™ comparable product code

Materials Chart






Sharpoint™ PLUS	Ethicon™	Syneture™	Material
PLAIN GUT	Surgical Gut	Plain Gut	Natural Absorbable Gut
CHROMIC GUT	Surgical Gut	Chromic Gut	Natural Absorbable Gut
POLYSYN FA™	Vicryl Rapide™	Velosorb™ Fast	Synthetic Absorbable, Braided, Coated Polyglycolic Acid
MONODERM™	Monocryl™	Caprosyn™/Biosyn™	Synthetic Absorbable, Monofilament Glycolide and E-caprolactone
POLYSYN™	Vicryl™	Polysorb™/Dexon™	Synthetic Absorbable, Braided, Coated Polyglycolic Acid
PDO	PDS™ II	Maxon™	Synthetic Absorbable, Monofilament Polydioxanone
SILK	Permahand Silk™	SofSilk™	Non-Absorbable, Braided, Coated Silk
NYLON™	Ethilon™	Monosof™/Dermalon™	Non-Absorbable, Monofilament Nylon
POLYPROPYLENE	Prolene™	Surgipro™II	Non-Absorbable, Monofilament Polypropylene
UNCOATED POLYVIOLENE™	Mersilene™	Surgidac™	Non-Absorbable, Braided, Uncoated Polyester



Sharpoint™ PLUS sutures are designed and manufactured to meet the exacting needs of surgeons worldwide.





A Guide To Surgical Sutures

Absorbable

Material	Material Description	BSR Profile	Mass Absorption
GUT  • Surgical Gut Suture - Plain/Chromic • Plain Gut/Chromic Gut	Twisted • Plain & Chromic • Packaged Dry or in Wetting Solution Common uses: General Closure, OB/Gyn, Urology, Bowel Anastomosis	In vivo strength retention: Varies	Varies depending upon type of suture (plain vs. chromic) and site of placement
POLYSYN FA™ (glycolic Acid)  • VICRYL RAPIDE™ • Velosorb™ Fast	Braided • Undyed (clear) • Coated Common uses: General Closure, OB/Gyn, Urology	In vivo strength retention: 48% at 7 days 0% at 14 days	Essentially complete between 49 and 63 days
MONODERM™ (glycolic Acid and e-caprolactone)  • MONOCRYL™ • Caprosyn™/Biosyn™	Monofilament • Undyed (clear) • Dyed (violet) Common uses: Skin Closure, OB/Gyn, Urology, Gastrointestinal	In vivo strength retention: 42% to 76% at 7 days 36% to 52% at 14 days	Essentially complete by 90 days
POLYSYN™ (glycolic Acid)  • VICRYL™ • Polysorb™/ Dexon™	Braided • Undyed (clear) • Dyed (violet) Monofilament • Undyed (clear) • Dyed (violet) • Coated & Uncoated Common uses: Skin Closure, OB/Gyn, Urology, Gastrointestinal	In vivo strength retention: 50% at 14 days 20% at 21 days	Essentially complete between 50 and 90 days
PDO (polydioxanone)  • PDS™ II • Maxon™	Monofilament • Dyed (violet) Common uses: Plastic Surgery, Orthopaedic, OB/Gyn, Gastrointestinal	In vivo strength retention: 80% to 90% at 14 days 60% to 82% at 28 days 47% to 79% at 42 days	Essentially complete between 180 and 220 days

• Comparable Ethicon™ Material * Comparable Syneture™ Material


Non-Absorbable

Material	Material Description	Qualities
SILK  • PERMA-HAND™ Silk * SofSilk™	Braided • Undyed (white) • Dyed (black) • Coated Common uses: General & Skin Closure, Gastrointestinal, Cardiovascular Surgery, Plastic Surgery, Ophthalmology	• Excellent handling and tying characteristics
NYLON™  • ETHILON™ • Monosof™/Dermalon™	Monofilament • Undyed (clear) • Dyed (blue) (black) • Packaged dry or in wetting solution Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Microsurgery, Cardiovascular Surgery, Gastrointestinal	• Passes easily through tissue
POLYPROPYLENE  • PROLENE™ • Surgipro™/Surgipro™II	Monofilament • Undyed (clear) • Dyed (blue) Common uses: General & Skin Closure, Neurosurgery, Plastic Surgery, Ophthalmology, Microsurgery, Cardiovascular Surgery	• Inert • Easy to handle • Passes easily through tissue
POLYVIOLENE™ (UNCOATED) (poly (ethylene-terephthalate))  • MERSILENE™ • Surgidac™	Braided • Undyed (white) • Dyed (green) • Uncoated Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Cardiovascular Surgery, Neurological Procedures	• Smooth • Easy to handle • Good knot security

• Comparable Ethicon™ Material * Comparable Syneture™ Material

SharpPoint™ PLUS sutures are available in a variety of absorbable and non-absorbable materials for use in a wide variety of surgical applications.

Needle Anatomy And Terminology



Needle Terminology

Length: Distance measured along the circumference of the needle, from the point to the swage. It is depicted on the packaging.

Shape: Different surgeries or surgical techniques demand different needle shapes.

Point: The sharpened part of the needle that first penetrates the tissue. The point consists of the area from the sharpened tip of the needle to the widest cross section of the needle body. The point determines ease of penetration and the initial size and shape of the hole.


Swage: Connection point of the suture and the needle. Micro teeth grip the suture material and prevent it from slipping out.

Wire Diameter: Gauge or thickness of the wire originally used to form the needle body.


Body: The part between the point and the swage of the needle.

Included Angle (Curvature): The fraction of the circle included in the curve of the needle. This is depicted on the packaging.


Needle Shapes And General Uses




Straight (aka: Keith): Skin




Bi-curve: Vasectomy Reversal, Ophthalmic




1/4 Circle: Ophthalmic



1/8 Circle: Ophthalmic



3/8 Circle: Skin Closure, Ophthalmic **1/2 Circle:** Skin Closure, GI, Cardiovascular, Ophthalmic



Needle Points And Cutting Edges

Needle Type	Point Geometry	Description	Commonly Used For
Conventional Cutting		Triangular point with cutting edge on the inner curvature	General skin closure, subcutaneous tissue, sometimes for ophthalmic surgery, plastic or reconstructive surgery
Precision Conventional Cutting		Triangular point with cutting edge on the inner curvature and precision-honed for added sharpness	General skin closure, subcutaneous tissue, sometimes for ophthalmic surgery, plastic or reconstructive surgery
Reverse Cutting		Triangular point with cutting edge on the outer curvature	Skin closure, ligamentous or fibrous tissues
Precision Reverse Cutting		Triangular point with cutting edge on the outer curvature and precision-honed for added sharpness	Primarily used in plastic surgery for delicate procedures and where a good cosmetic result is desired
Taper Point		Round shaft, straight or curved, taper point with no cutting edge	Soft tissue closure such as fascia, vascular, gastrointestinal, and most soft tissue below the skin
Taper Cutting		Four (4) cutting edges with a taper shaft	For penetration through tough tissue
Side Cutting Lancet		Spatula needle with the cutting edge on the inner, or outer curvature	Ophthalmic and microsurgery
Center Point Spatula		Two cutting edges in a horizontal plane	Ophthalmic surgery for muscle and retinal repair. Also used for delicate plastic surgery
Straight Cutting		Reverse cutting tip on a straight shaft	General skin closure
M.E.T.™ Point		Triangular faceted point tapering to a round body	Microsurgery
Diamond Point™		Modified taper cutting needle	Penetration of tougher tissues like the dermis
UltraGlide™		Six (6) faceted sides and two (2) cutting edges	Ophthalmic surgery for cornea and sclera closure
Vas Cutting		Cutting edges aid in penetration of multiple layers of tissue	Microsurgery and urology

Needle Chart

Center Point Spatula Needles

Surgical Specialties Ethicon™ Syneture™

1/4 Circle

DO5	RD-1	SS-24	
-----	------	-------	--

Conventional Cutting Needles

Surgical Specialties Ethicon™ Syneture™

3/8 Circle


DSC80	GCC-90	
-------	--------	--

Diamond Point Needles

Surgical Specialties Ethicon™ Syneture™

3/8 Circle

DE7	G-6	HE-6	
-----	-----	------	--



M.E.T.™ Point Needles

Surgical Specialties Ethicon™ Syneture™

Bi-curve

BRM5			
VRM4	V75-3		
VRM5			
VRM6			
DRM4		MVK-70-3	
DRM5		MVK-100-4	
DRM6			
HRM4			

Precision Conventional Cutting Needles

Surgical Specialties Ethicon™ Syneture™

3/8 Circle

DSCM13	PC-1	P-13	
DSCM16	PC-3	PC-11	
DSCM19	PC-5	P-12	

Needle Chart

Precision Reverse Cutting Needles

Surgical Specialties	Ethicon™	Syneture™
3/8 Circle		
DSM7	G-6	HE-6
DSM11	P-1	P-10
DSM13	P-3	P-13
DSM16	PS-3	P-11
DSM18	PS-2	P-12
DSM19	PS-2	P-12
DSM24	PS-1	P-14
1/2 Circle		
HSM8	G-7	HE-7

Side Cutting Lancet Needles

Surgical Specialties	Ethicon™	Syneture™
Bi-curve		
BSL5	TG6-C	•
1/8 Circle		
SSL5	•	•
SSL15	•	•
1/4 Circle		
VSL7	TG100-8	SE100-8
VSP8	S-29, S-24	SD-1, SS-24
VCL8	S-14	SS-14
3/8 Circle		
DSL6	TG140-6	SE140-6
DSL7	TG140-8	SE140-8
1/2 Circle		
HSL5	TG160-6	SE160-6
HSL6	TG160-8	SE160-8
HSL7	•	•

Reverse Cutting Needles

Surgical Specialties	Ethicon™	Syneture™
3/8 Circle		
•	C-2	C-1
DS12	C-2	C-1
DS14	C-3	C-12
DS16	FS-3	C-12
DS18	FS-2	C-13
DS24	FS-1	C-14
DS26	FS	C-15
DS30	FSL	C-16
DS36	FSLX	C-17
DS40	FSLX	C-17
DS90	•	•
1/2 Circle		
•	M-1	C-21
HS16	•	•
HS18	•	•
•	J-1	C-22
HS24	X-1	C-23
•	CP-2	GS-10
•	CP-1	GS-11
•	CP	GS-12

Needle Chart

Taper Cutting Needles

Surgical Specialties	Ethicon™	Syneture™
3/8 Circle		
DRM4	BV-100-3	MVK-70-3
DRM5	BV-100-4	•
DRM6	BV-130-5	MV-135-5
•	CC-1	•
•	RB-1	CV-23
1/2 Circle		
•	V-5	KV-5
HRT17	V-5	KV-5
HRT31	V-34	KV-34
HRT36	V-34	KV-34
HRT38	V-34	KV-34

Taper Point Needles

Surgical Specialties	Ethicon™	Syneture™
1/4 Circle		
VR5	V75-3	•
3/8 Circle		
DR4	BV75-3	MV70-3
DR5	BV130-4	MV135-4
DR6	BV130-5	MV135-5
DR12	C-1	CV-1
DR20	SH-1	CV-24
1/2 Circle		
HR13	TF	CV-22
HR18	RB-1	CV-23
HR22	SH-1	CV-25
HR26	SH	V-20
•	CT-2	GS-22
HR30	•	V-30
HR36	MH	V-26
HRX36	CT-1	GS-21
HRX40	CT	GS-24X

Straight Cutting Needles

Surgical Specialties	Ethicon™	Syneture™
GS60	KS	SC-2

UltraGlide™ Needles



Surgical Specialties	Ethicon™	Syneture™
Bi-curve		
BGL5	CSB-6	•
1/4 Circle		
VGL7	CS90-6	•
3/8 Circle		
DGL6	CS140-6	SE-140-6
DGL7	CS140-6	SE-140-6
1/2 Circle		
HGL5	CS160-6	SE-160-6
HGL7	CS175-6	•

Vas Cutting Needles






Surgical Specialties	Ethicon™	Syneture™
1/2 Circle		
HSV6	•	•

Plain Gut

Precision Conventional Cutting

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSCM11 11mm 3/8 Circle  PC-1 P-13	Plain Gut	18"/45cm	B1415N •						
DSCM13 13mm 3/8 Circle  PC-1 P-13	Plain Gut	18"/45cm	B2711N •						

Precision Reverse Cutting

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSM11 11mm 3/8 Circle  P-1 P-10	Plain Gut	12"/30cm	B2698N •						
		18"/45cm	B660N/ B774N 774G SG1645						
DSM11 11mm 3/8 Circle  P-1 P-10	Plain Gut	18"/45cm	B770N 770G SG1645						
DSM12 12mm 3/8 Circle  C-2 C-1	Plain Gut	10"/45cm		B521NS •	B531NS •				
		18"/45cm		B686N 686G SG5686	B536N • SG721				
DSM13 13mm 3/8 Circle  P-3 P-13	Plain Gut	18"/45cm			B1644N 1644G SG5654G				
DS18 19mm 3/8 Circle  FS-2 C-13	Plain Gut	18"/45cm		B0500N • SG751	B2463N • SG752	B592N • SG5630G			
		27"/70cm		H820G SG820G	B821N H821H SG821G	B822N H822H SG822G			
DS24 24mm 3/8 Circle  FS-1 C-14	Plain Gut	27"/70cm				B810N H810H SG810			

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Plain Gut

Reverse Cutting







Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DS24 24mm 3/8 Circle  FS-1 C-14	Plain Gut	27"/70cm						B823N H823H SG823	
C26 15mm 1/2 Circle  M-1 C-21	Plain Gut	18"/45cm			B537N •				
HSM22 22mm 1/2 Circle  X-1 C-23	Plain Gut	18"/45cm					B612N 612G SG612G		
C21 26mm 1/2 Circle  CP-2 GS-10	Plain Gut	18"/45cm				B0599N • SG1633G	B0593N • SG612G		



Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Chromic Gut







Precision Reverse Cutting ▾

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
C1 P-1 P-10 11mm 3/8 Circle 	Chromic Gut	12"/30cm	C796N 796G G1780						
C3 P-3 P-13 13mm 3/8 Circle 	Chromic Gut	18"/45cm		C687N 687G SG5687G	C1654N 1654G SG5644G				
C3 P-3 DX-13 13mm 3/8 Circle 	Chromic Gut	18"/45cm		C792N 792G •					
DSM16 PS-3 P-11 16mm 3/8 Circle 	Chromic Gut	18"/45cm		C1636N 1636G •					
DSM19 PS-2 P-12 19mm 3/8 Circle 	Chromic Gut	18"/45cm			C558N 1637G SG5637G/ SG675	C525N • SG676			
		27"/70cm		C634N 634G •	C635N 635H SG635G	C1638N/ C636N 1638H/ 636H SG5163G/ SG636G			
		30"/75cm		C0935N • SG634G					
DSM24 FS-1 C-14 24mm 3/8 Circle 	Chromic Gut	27"/70cm						C637N 637H SG637G	

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Chromic Gut






Reverse Cutting ▾

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
C16 C-2 C-1 11mm 3/8 Circle 	Chromic Gut	12"/30cm		C543N • •					
C17 C-2 C-1 12mm 3/8 Circle 	Chromic Gut	12"/30cm			752G SG5644G				
DS24 FS-1 C-14 24mm 3/8 Circle 	Chromic Gut	30"/75cm				C100N • •			
DS26 FS C-15 26mm 3/8 Circle 	Chromic Gut	27"/70cm						C0562N • •	
HS15 M-1 C-21 15mm 1/2 Circle 	Chromic Gut	18"/45cm		C1642N 1642G •	C1653N 1653G •				
HS20 J-1 C-22 20mm 1/2 Circle 	Chromic Gut	18"/45cm			C1643N 1643G •				

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Chromic Gut


Taper Point

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
T28 18mm 1/2 Circle  RB-1 CV-23	Chromic Gut	14"/35cm		C0590N • UG202					
HR22/ T4 22mm 1/2 Circle  SH-1 CV-25	Chromic Gut	14"/35cm			C181N G181H GG181				
		30"/75cm				C182N G182H GG182			
HR26 26mm 1/2 Circle  SH V-20	Chromic Gut	30"/75cm			C121N G121H GG121	C122N 6122H GG122	C123N 6123H GG123		
HR30 30mm 1/2 Circle  V-20	Chromic Gut	30"/75cm					C0244N • GG130		
HR36 36mm 1/2 Circle  MH V-26	Chromic Gut	30"/75cm						C128N G128H GG128	


Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

PolySyn FA™ / Polyglycolic Acid


Precision Reverse Cutting

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSM11 11mm 3/8 Circle  P-1 P-10	Undyed Braided	18"/45cm		M490N VCP490G SL1689G					
DSM13 13mm 3/8 Circle  P-3 P-13	Undyed Braided	18"/45cm		M493N VR493 •	M494N VR494 SL5691G				
		19mm 3/8 Circle  PS-2 P-12	18"/45cm			M496N VR496 •	M497N VR497 •		
DS24 24mm 3/8 Circle  FS-1 C-14	Undyed Braided	27"/70cm			M426N VR426 •	M0423N • •			
		24mm 1/2 Circle  X-1 C-23	18"/45cm				M442N • VC442H		
HS24 24mm 1/2 Circle  X-1 C-23	Undyed Braided	18"/45cm					M456N • CL53MG		

Taper Cutting




Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
HRT36 36mm 1/2 Circle  V-34 KV-34	Undyed Braided	36"/90cm						M917N VR917 J9962T	

Taper Point




Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
HR17 17mm 1/2 Circle  RB-1 CV-23	Undyed Braided	27"/70cm			M214N VR214 •				

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Precision Conventional Cutting ▲

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSCM13 13mm 3/8 Circle PC-1 P-13 	Undyed Monofilament	18"/45cm	Y833N •						
DSCM16 16mm 3/8 Circle PC-3 PC-11 	Undyed Monofilament	10"/25cm		Y844NS Y844G SC5626G	Y845NS Y845G SC5686G				
		18"/45cm		Y844N Y844G SC5626G	Y845N Y845G SC5686G				
DSCM19 19mm 3/8 Circle PC-5 P-12 	Undyed Monofilament	10"/25cm			Y823NS •				
		18"/45cm		Y822N Y822G SC5616G	Y823N Y823G SC5617G				

Precision Reverse Cutting ▼



Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSM11 11mm 3/8 Circle P-1 P-10 	Undyed Monofilament	10"/25cm		Y490NS •					
		18"/45cm	Y489N Y489G SC3625G	Y490N Y490G SC3626G					
DSM13 13mm 3/8 Circle P-3 P-13 	Undyed Monofilament	10"/25cm		Y493NS •	Y494NS •				
		18"/45cm	Y492N Y492G SC5688G	Y493N Y493G SC5689G/SC5687	Y494N Y494G SC5690G				
	Violet Monofilament	18"/45cm		Y463N Y463G SC5587G	Y464N Y464G SC5590G				
DSM19 19mm 3/8 Circle PS-2 P-12 	Undyed Monofilament	10"/25cm		Y495NS •	Y426NS •	Y427NS •			
		18"/45cm		Y495N Y495G SC5616G	Y496N Y496G SC5635G	Y497N Y497G SC5618G			
		27"/70cm			Y426N Y426G SC5637/SC5637G	Y427N Y427H •			
	Violet Monofilament	18"/45cm			Y513N Y513G SC5627G				

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Precision Reverse Cutting (cont'd) ▼

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSM24 24mm 3/8 Circle PS-1 P-14 	Undyed Monofilament	18"/45cm			Y682N Y682H SC5679	Y936N Y936H SC5780G			
HSM9 9mm 1/2 Circle P-2 P-21 	Undyed Monofilament	18"/45cm		Y503N •	Y504N •				

Taper Point ⊙

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
HR13 13mm 1/2 Circle TF CV-22 	Undyed Monofilament	27"/70cm	Y432N Y432H UM886	Y433N Y433H SM432	Y434N Y434H SM431				
HR26 26mm 1/2 Circle CT-2 NCT-2 	Undyed Monofilament	27"/70cm						Y417D •	



Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

PolySyn™/Polyglycolic Acid




Diamond Point

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DE7 8.13mm 3/8 Circle  G-6 HE-6	Violet Braided	12"/30cm	G551N J551G L1551K						

Precision Conventional Cutting

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSCM13 13mm 3/8 Circle  PC-1 P-13	Undyed Braided	18"/45cm		G834N J834G SL5687G					
DSCM16 16mm 3/8 Circle  PC-3 PC-11	Undyed Braided	18"/45cm			G835N J845G SL1657				





Precision Reverse Cutting

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSM11 11mm 3/8 Circle  P-1 DX-11	Undyed Braided	8"/20cm	G489NS •						
		18"/45cm		G490N J490G SL1698G					
DSM13 13mm 3/8 Circle  P-3 P-13	Undyed Braided	8"/20cm	G492NS •						
		10"/25cm		G493NS •	G494NS •				
		18"/45cm	G492N J492 SL5688G	G493N J493G SL5687G	G494N J494G SL5690 / SL5691G				
	Violet Braided	18"/45cm		G463N J463G SL5687G	G464N J464G SL5691G				
DSM16 16mm 3/8 Circle  PS-3 P-11	Undyed Braided	18"/45cm	G499N J499G XL843	G500N J500G SL5626G			G0508N • SL5621		

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

PolySyn™/Polyglycolic Acid




Precision Reverse Cutting (cont'd)

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSM19 19mm 3/8 Circle  PS-2 P-12	Undyed Braided	8"/20cm		G495NS •	G496NS •				
		18"/45cm	G794N •	G495N J495G XL822/ SL5636G	G496N J496G XL823/ SL5627G	G497N XL824/ SL5628G			
		27"/70cm		G421N J421H SL690G	G426N/ G422N/ G422D J426H/ J422H SL5637/ SL691G	G427N/ G423N J427H/ J423H SL5638/ SL822G			
		30"/75cm			G691D •				
	Violet Braided	18"/45cm				G393N J393H SL822			
		27"/70cm			G0392N J392H •				
DSM24 24mm 3/8 Circle  PS-1 P-14	Undyed Braided	18"/45cm					G683N J683G SL5640G		
HSM9 9mm 1/2 Circle  P-2 P-21	Undyed Braided	18"/45cm		G503N J503G SL1644G	G504N J504G SL1645G				
HSM16 16mm 1/2 Circle  PS-4 P-12	Undyed Braided	18"/45cm		G506N J506G SL5632G	G507N J507G SL5633G				


Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

PolySyn™/Polyglycolic Acid

Reverse Cutting ▾

Needle	Material	Length	Size					
			5-0	4-0	3-0	2-0	0	1
DS18 19mm 3/8 Circle  FS-2 C-13	Violet Braided	27"/70cm	G391N J391H SL634	G392N J392H SL635G				
DS24 24mm 3/8 Circle  FS-1 C-14	Undyed Braided	27"/70cm		G441N J441G SL692G	G442N J442H SL693	G443N/ G443D J443H SL823		
	Violet Braided	27"/70cm		G451N • •				
HS40 40mm 1/2 Circle  CP GS-12	Violet Braided	27"/70cm						G480N J480H CL805
		36"/90cm						J486H •

Side Cutting Lancet ▽

Needle	Material	Length	Size			
			6-0	5-0	4-0	3-0
VSL7 1/4 Circle Length 6.60mm Curvature 94  TG100-8 SE-100-8	Violet Braided	18"/45cm	G544N J544G L1750K			
VCL8 1/4 Circle Length 8.51mm Curvature 112  SG-14 SS-14	Undyed Braided	18"/45cm	G670N J670G L1670	G671N J671G L1669		
	Violet Braided	8"/20cm	G590N J590G L1740K			



Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

PolySyn™/Polyglycolic Acid

Side Cutting Lancet (cont'd) ▽

Needle	Material	Length	Size			
			6-0	5-0	4-0	3-0
VSP8 1/4 Circle Length 7.62mm Curvature 110  S-29 SD-1	Violet Braided	12"/30cm	G556N J556G L1446K			
		18"/45cm	G555N J555G •			
VSP08 1/4 Circle Length 8.00mm Curvature 90  S-24 SS-24	Violet Braided	12"/30cm	G552N J552G L1742K	G553N J553G L1785K		







Taper Cutting ▽

Needle	Material	Length	Size					
			6-0	5-0	4-0	3-0	2-0	0
HRT17 17mm 1/2 Circle  V-5 KV-5	Violet Braided	18"/45cm		G0502N • •				
HRT31 31mm 1/2 Circle  V-5 KV-5	Violet Braided	30"/75cm						G9391N W9391 •

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

PolySyn™/Polyglycolic Acid




Taper Point

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
HR13 13mm 1/2 Circle  TF CV-22	Undyed Braided	27"/70cm		G433N •					
	Violet Braided	18"/45cm	G981N •						
HR17 17mm 1/2 Circle  RB-1 CV-23	Undyed Braided	27"/70cm			G714N J714D UL212				
	Violet Braided	30"/75cm			G0304N •				
T14 22mm 1/2 Circle  SH-1 CV-24	Violet Braided	27"/70cm			G310N J310H GL181	G311N J311H GL182			
HR26 26mm 1/2 Circle  SH V-20	Undyed Braided	27"/70cm			G416N J232H GL322/ GLS322				
	Violet Braided	27"/70cm			G316N J316H GL122	G317N J317H GL123			
T5 26.5mm 1/2 Circle  CT-2 GS-22	Undyed Braided	27"/70cm			G232N J232H CL829	G269N J269H CL830			
	Violet Braided	27"/70cm			G332N J332H CL882	G333N J333H CL883			
HR30 30mm 1/2 Circle  • V-30	Violet Braided	30"/75cm			G0322N •	G323N •			


Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

PolySyn™/Polyglycolic Acid

Taper Point (cont'd)


Needle	Material	Length	Size						
			5-0	4-0	3-0	2-0	0	1	
T8 36.6mm 1/2 Circle  CT-1 GS-21	Undyed Braided	18"/45cm				G839N J839D CL13MG			
		27"/70cm				G260N J260H CL844			
		36"/90cm				G945N J945H CL953			
T10 40mm 1/2 Circle  CT GS-24	Violet Braided	27"/70cm				G339N J339H CL930	G340N J340H CL812	G341N J341H CL813	
		36"/90cm				G354N J345H •	G0518N •		
T10 40mm 1/2 Circle  CT GS-24	Violet Braided	27"/70cm					G352N J352H CL802	G353N J353H CL803	
		36"/90cm					G358N J358H •	G359N J359H •	

Reverse Cutting / Taper Point




Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
HSM36/ 27mm/36.6mm HR36 1/2 Circle  CP-1 GS-11 CT-1 GS-21	Violet Braided	54"/135cm						J532H •	

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".








Precision Reverse Cutting ▾

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSM19 19mm 3/8 Circle  PS-2 P-12	Violet Monofilament	18"/45cm			R513N Z513G				



Taper Point ⊙

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
HR17 17mm 1/2 Circle  RB-1 NRB-1	Violet Monofilament	27"/70cm				R305D Z305H			
T5 26mm 1/2 Circle  SH/CT-2 NSH/NC-2	Violet Monofilament	27"/70cm				R316D Z316H	R317D/ R333D Z317H/ Z333H	R334D Z334H	
HR3 36mm 1/2 Circle  CT-1 NCT-1	Violet Monofilament	27"/70cm					R339D Z339H	R340D Z340H	

Precision Reverse Cutting ▾

Needle	Material	Length	Size						
			7-0	6-0	5-0	4-0	3-0	2-0	
DSM7 7mm 3/8 Circle  G-6 HE-6	Black Braided	12"/30cm	• SL1768K						
DSM7 7mm 3/8 Circle  G-6 HE-6	Black Braided	18"/45cm	D0191N 768G SL1768K	D769N 769G					
DSM11 11mm 3/8 Circle  P-1 P-10	Black Braided	18"/45cm		D639N 639G SS1639G					
DSM11 11mm 3/8 Circle  P-1 P-10	Black Braided	18"/45cm		780G S1750					
C3 13mm 3/8 Circle  P-3 P-13	Black Braided	18"/45cm			D640N 640G SS5640G	D641N 641G SS5641G			
PC31 19mm 3/8 Circle  PS-2 P-12	Black Braided	18"/45cm				D1677N 1677G SS5677G			
HSM8 8mm 1/2 Circle  G-7 HE-7	Black Braided	18"/45cm		D765N 765G S1675					








Reverse Cutting ▾

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DS12 12mm 3/8 Circle  C-2 C-1	Black Braided	18"/45cm	D1639N 1639G S1172	D0374N • S1173	D736N 736G S1174				
DS16 16mm 3/8 Circle  FS-3 C-12	Black Braided	18"/45cm	D0750N • SS681	D1678N 1678H SS5678	D0239N • SS686	D0389N •			
		30"/75cm			D0208N •				

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".


Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Reverse Cutting (cont'd) ▾



Needle	Material	Length	Size					
			5-0	4-0	3-0	2-0	0	1
DS18 19mm 3/8 Circle  FS-2 C-13	Black Braided	18"/45cm	D682N 682G SS682G	D683N 683G SS683G	D1679N 1679H SS694/ SS5679			
		30"/75cm			D0768N • SS745			
DS24 24mm 3/8 Circle  FS-1 C-14	Black Braided	18"/45cm		D629N 629G SS629	D684N 684G SS684G	D0238N 685G SS695		
		30"/75cm			D0775N • •			
DS26 26mm 3/8 Circle  FS C-15	Black Braided	18"/45cm				D685N 685G SS685G		
		30"/75cm			D0212N • •			
DS30 30mm 3/8 Circle  FSL C-16	Black Braided	18"/45cm			D0765N • •		D678N 678G •	
		30"/75cm			D675N 675H SS675	D679N 679H SS679	D680N 680H SS680	
DS36 36mm 3/8 Circle  FSLX C-17	Black Braided	30"/75cm				D0213N • SS648		D0423N • SS646
C26 15mm 1/2 Circle  M-1 C-21	Black Braided	18"/45cm	D779N • •	D777N • •	D778N • •			
HS16 16mm 1/2 Circle  ••	Black Braided	18"/45cm			D2460N • •			

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Reverse Cutting (cont'd) ▾

Needle	Material	Length	Size				
			6-0	5-0	4-0	3-0	2-0
HSM17 17mm 1/2 Circle  J-1 C-22	Black Braided	18"/45cm			734G SS733		
HS18 18mm 1/2 Circle  ••	Black Braided	18"/45cm				D0783N • •	
HSM22 22mm 1/2 Circle  ••	Black Braided	18"/45cm				D632N 632G SS732	
HS26 26mm 1/2 Circle  CP-2 GS-10	Black Braided	18"/45cm				D782N • •	


Side Cutting Lancet ▾

Needle	Material	Length	Size		
			8-0	7-0	6-0
VSP8 1/4 Circle Length 8mm Curvature 90  TG100-8 SE-100-8 Radius 5.08mm Chord 7.24mm Wire .35mm	Black Braided	18"/45cm			D1780N 1780G S2780
DSL7 3/8 Circle Length 6.15mm Curvature 135  TG140-8 SE-140-8 Radius 2.54mm Chord 4.83mm Wire .15mm	Black Braided	12"/30cm	D2526N 599G •		
		18"/45cm		D7733N 7733G S1768K	


Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
 •• indicates that other manufacturers do not have a comparable product code or needle.
 N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Silk




Straight Cutting ▾

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
GS60 60mm Straight  KS SC-2	Black Braided	30"/75cm						D623N 623H SS623	

Taper Cutting Ⓞ

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DRT14 13mm 3/8 Circle  CC-1 KZ-11	Black Braided	18"/45cm	D7729N 7729G						

Taper Point Ⓞ

Needle	Material	Length	Size							
			6-0	5-0	4-0	3-0	2-0	0	1	
DR12 13mm 3/8 Circle  C-1 CV-1	Black Braided	18"/45cm	D707N 707G VS802		D7734N 7734G					
HR26 26mm 1/2 Circle  SH V-20	Black Braided	30"/75cm				D832N K832H GS832				
HR36 36mm 1/2 Circle  MH V-26	White Braided	30"/75cm								D2543N CS745


Pre-Cut Sutures

Needle	Length	Size				
		0	1	2	3	4
Silk, Black Braided Strands	8 - 24"/60cm					D2578N • •
	13 - 24"/60cm	D76N • •		D0492N • •		

Product numbers are color coded respectively as Sharpoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Nylon




Conventional Cutting △

Needle	Material	Length	Size					
			6-0	5-0	4-0	3-0	2-0	
DSC80 80mm 3/8 Circle  GCC-90	Black Monofilament Pliable	40"/100cm						A2555N • CN791

Precision Conventional Cutting △




Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
DSCM16 6mm 3/8 Circle  PC-3 PC-11	Black Monofilament	18"/45cm	1866G SN1966	A1865N 1865G SN3965	A1864N 1864G SN1964				
	Black Monofilament Pliable	18"/45cm		A1965N 1965G SN3965					

Precision Reverse Cutting ▽



Needle	Material	Length	Size				
			7-0	6-0	5-0	4-0	3-0
DSM11 11mm 3/8 Circle  P-1 P-10	Black Monofilament	8"/20cm		• •			
		18"/45cm	A1696N 1696G SN3686	A697N 697G SN3697			
DSM12 12mm 3/8 Circle  C-2 C-1	Black Monofilament	10"/25cm				A931NS • •	
DSM13 13mm 3/8 Circle  P-3 P-13	Black Monofilament	8"/20cm		A1273NS/A1698NS • •	A1284NS • •	A699NS • •	
		18"/45cm		A1698N 1698G SN5696G	A698N 698G SN5698G	A699N 699G SN5699G	
	Clear Monofilament	18"/45cm			A690N 690G SN5690G	691G SN5691G	

Product numbers are color coded respectively as Sharpoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Precision Reverse Cutting (cont'd) ▾








Needle	Material	Length	Size				
			7-0	6-0	5-0	4-0	3-0
DSM16 16mm 3/8 Circle  PS-3 P-11	Black Monofilament	18"/45cm		1665G SN5767G	A1668N 1658G SN5668G		
DSM18 18mm 3/8 Circle  PS-2 P-12	Black Monofilament	8"/20cm			A1666NS •	A1667NS •	
	Black Monofilament	18"/45cm		A1281N • SN5665G	A1666N 1666G SN5666G	A1667N 1667G SN5667G	A1669N/ A1669D 1669G SN5669G
DSM24 24mm 3/8 Circle  PS-1 P-14	Clear Monofilament	18"/45cm				A1611N 1611G SN5677	A2506N • SMM5044
	Black Monofilament	18"/45cm					A1663N 1663G SN5633G
	Black Monofilament	30"/75cm					A669N/ A669D 669H SN669G

Reverse Cutting ▾

Needle	Material	Length	Size				
			6-0	5-0	4-0	3-0	2-0
DS12 12mm 3/8 Circle  C-2 C-1	Black Monofilament	9"/23cm	A2557NS •	A2558NS •			
		10"/23cm	A911NS •	A913NS •			
		18"/45cm	A667N 667G SN667G	A668N 668G SN871	A921N 1954G SN872		
	Clear Monofilament	10"/25cm		A668NS •			
DS16 16mm 3/8 Circle  FS-3 C-12	Black Monofilament	18"/45cm	A0181N 1655G SN660G	A0556N 1668G SN659G			


Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Reverse Cutting (cont'd) ▾


Needle	Material	Length	Size						
			5-0	4-0	3-0	2-0	0	1	
DSM19 19mm 3/8 Circle  FS-2 C-13	Black Monofilament	10"/25cm	A661NS •	A662NS •					
		18"/45cm	A661N/ A661D/ 661G/ 661H SN661G	A662N/ A662D 662G SN662G					
DS24 24mm 3/8 Circle  FS-1 C-14	Black Monofilament	10"/25cm		A947NS 1629H SN658					
		18"/45cm		A1629N/ A1629D 1629H SN658	A663N 663G SN663G				
DS26 26mm 3/8 Circle  FS C-15	Black Monofilament	18"/45cm				A664N 664G SN764G			
	Blue Monofilament	30"/75cm		A1328N •	A8685N 8685H SP868G				
C10 30mm 3/8 Circle  FSL C-16	Clear Monofilament	27"/70cm		A1671N 1671H SN671					
DS40 40mm 3/8 Circle  FSLX C-17	Blue Monofilament	30"/75cm				A0616N •	A0620N •		
DS90 90mm 3/8 Circle  ••	Blue Monofilament	40"/100cm					A0637N •		
HSM16 16mm 1/2 Circle  PC-12 DHX-16	Black Monofilament	18"/45cm		A1662N 1662G SN5662G					

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Straight Cutting ▽




Needle	Material	Length	Size				
			6-0	5-0	4-0	3-0	2-0
SC2 60mm Straight  KS SC-2	Black Monofilament	30"/75cm				A627N 627H SN627	

Straight Taper Point ⊙






Needle	Material	Length	Size				
			6-0	5-0	4-0	3-0	2-0
S 61mm Straight  KS SC-2	Black Monofilament	30"/75cm					A933N •

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Precision Conventional Cutting ▲

Needle	Material	Length	Size					
			6-0	5-0	4-0	3-0	2-0	0
CP1 11mm 3/8 Circle  PC-1 P-13	Blue Monofilament	8"/20cm	J1057NS •	J1058NS •				
DSCM13 13mm 3/8 Circle  PC-1 P-13	Blue Monofilament	18"/45cm		J8618N •				
DSCM16 16mm 3/8 Circle  PC-3 P-11	Blue Monofilament	18"/45cm		J2451N/ J8635N 8635G	J8634N 8634G			





Precision Reverse Cutting ▽

Needle	Material	Length	Size				
			7-0	6-0	5-0	4-0	3-0
DSM11 11mm 3/8 Circle  P-1 P-10	Blue Monofilament	8"/20cm		J8697NS •			
		18"/45cm	J8696N 8696G SP1696G	J8697N 697G SN3697			
DSM13 13mm 3/8 Circle  P-3 P-13	Blue Monofilament	8"/20cm		J1077NS •	J1078NS •	•	
		18"/45cm		J8695N 8695G/ 8617G SP5695G/ SP1617	J8698N 8698G SP5698G	J8699N 8699G SP5699G	
DSM16 16mm 3/8 Circle  P-3 P-13	Blue Monofilament	18"/45cm		J8680N 8680G SP5680G	J8681N 8681G SP5681G		
DSM18 18mm 3/8 Circle  PS-2 P-12	Blue Monofilament	8"/20cm			J8686NS •	J8682NS •	
		18"/45cm			J8686N 8686G SP5686G	J8682N 8682G SP5682G	J8687N 8687G SP5687G
DSM19 19mm 3/8 Circle  PS-2 P-12	Blue Monofilament	18"/45cm			J661N 8661G SP661G	J683N • SP683G	
DSM24 24mm 3/8 Circle  PS-1 P-14	Blue Monofilament	18"/45cm					J8663N 8663G SP5663G


Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Polypropylene


Reverse Cutting ▾

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
C17 12mm 3/8 Circle  C-2 C-1	Blue Monofilament	10"/25cm	J1011NS •	J1013NS •	J1031NS •				
DS14 14mm 3/8 Circle  C-3 C-12	Blue Monofilament	18"/45cm	J1712N • SP660G						
DS18 18mm 3/8 Circle  FS-2 C-13	Blue Monofilament	18"/45cm		J8661N 8661G SP661G	J8683N J8683D 8683G SP683G	J8665N 8665G SP665G			
DS24 24mm 3/8 Circle  FS-1 C-14	Blue Monofilament	18"/45cm			J8629N 8629H •	J8684N 8684G SP684G			
DS26 26mm 3/8 Circle  FS C-15	Blue Monofilament	18"/45cm				J1239N • SP684G	J8685N 8685G SP868G		

Straight Cutting ▾

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
GS60 60mm Straight  KS SC-2	Blue Monofilament	30"/75cm						J8623N 8623H SP623	

Taper Point ⊙

Needle	Material	Length	Size						
			8-0	7-0	6-0	5-0	4-0	3-0	0
HRX40 40mm 1/2 Circle  MH V-26	Blue Monofilament	5"/13cm							J8424N 8424H •


Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Polyviolene™ / Polyester - Uncoated


Reverse Cutting ▾

Needle	Material	Length	Size						
			7-0	6-0	5-0	4-0	3-0	2-0	
C22 16mm 3/8 Circle  FS-3 C-12	Green Braided Uncoated	18"/45cm			F681N X681G •				
DS18 18mm 3/8 Circle  FS-2 C-13	Green Braided Uncoated	18"/45cm				F683N X683G •			
	White Braided Uncoated	18"/45cm				F633N R633G XT692			
DS24 24mm 3/8 Circle  FS-1 C-14	White Braided Uncoated	18"/45cm						F647D R647H •	
DS26 26mm 3/8 Circle  FS C-15	White Braided Uncoated	18"/45cm							F665D R665H •

Side Cutting Lancet ◀▶

Needle	Material	Length	Size						
			7-0	6-0	5-0	4-0	3-0	2-0	
VSP8 8mm 1/4 Circle  S-24 SS-24	White Braided Uncoated	18"/45cm			F893N W893 •				

Taper Cutting ⊖

Needle	Material	Length	Size						
			6-0	5-0	4-0	3-0	2-0	0	
HRT17 17mm 1/2 Circle  V-5 KV-5	Green Braided Uncoated	12"/30cm Loop		F951N X951H •	F2600N •	F872N •			

Product numbers are color coded respectively as SharpPoint™ PLUS (dark blue), comparable Ethicon™ (black), and comparable Syneture™ (green).
•• indicates that other manufacturers do not have a comparable product code or needle.
N=box of 12; D=box of 36; S=SmallStitch™ suture lengths of 8"-10".

Visit SurgicalSpecialties.com and Click “Find A Product” for SKUs and Cross-Reference Information.

Surgical Specialties Corporation, 50 Braintree Hill Office Park, Suite 101, Braintree, MA 02184, U.S.A.
Phone: 610.404.3520 • 877.991.1110 (USA and Canada only) • Fax: 610.404.4010
All trademarks are property of their respective owners. ©2016 Surgical Specialties Corporation.
All Rights Reserved. SSC-028 R2 4/16

**SURGICAL
SPECIALTIES**
CORPORATION™