

Products designed with a purpose...
...to help you help them.



SURFLO® Winged Infusion Sets

Take a look at our wings versus the competition and you'll see why SURFLO is a very well recognized brand. Terumo's soft, flexible interlocking wings readily conform to the body's contour. The wing is angled for ease of positioning of the needle during venipuncture, for less risk of vessel damage or occlusion of the lumen.


The precisely honed needle with ultra-thin wall design assures ease of penetration and continued access with a single insertion.

LATEX FREE

 **TERUMO®**

TERUMO MEDICAL PRODUCTS

Surflo® Winged Infusion Sets



Intermittent Plug, 3 1/2" Tubing

Packed: 100/box, 5 boxes/case (500)

Product Code	Description
SV*19BI	19G X 3/4" Thin Wall Needle,
SV*21BI	21G x 3/4" Ultra Thin Wall Needle
SV*23BI	23G x 3/4" Ultra Thin Wall Needle

8" Tubing

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
SV*25EL	25G x 1/2" Regular Wall Needle
SV*27EL	27G x 1/2" Regular Wall Needle

12" Tubing

Packed: 100/box, 10 boxes/case (1,000)

Product Code	Description
SV-18BLK	18G x 3/4" Thin Wall Needle
SV-19BLK	19G x 3/4" Thin Wall Needle
SV-21BLK	21G x 3/4" Ultra Thin Wall Needle
SV-22BLK	22G x 3/4" Ultra Thin Wall Needle
SV-23BLK	23G x 3/4" Ultra Thin Wall Needle
SV-25BLK	25G x 3/4" Thin Wall Needle

3 1/2" Tubing

Packed: 100/box, 5 boxes case (500)

Product Code	Description
SV*19BLS	19G x 3/4" Thin Wall Needle
SV*21BLS	21G x 3/4" Ultra Thin Wall Needle
SV*23BLS	23G x 3/4" Ultra Thin Wall Needle
SV*25BLS	25G x 3/4" Thin Wall Needle


TerumoTMP.com/Infusion

US Office

2101 Cottontail Lane
Somerset, NJ 08873
Phone: 732.302.4900
Fax: 732.302.3093
Toll free: 800.283.7866
Fax for Orders: 800.411.5870

Canada Office

Terumo Medical Corporation
C/O UPS - Supply Chain Solutions
1165 Creditstone Road
Vaughan, Ontario L4K 4N7
Phone: 800.668.5489
Fax: 905.760.1215

 and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. Surflo is a trademark owned by Terumo Medical Corporation, Elkton, MD, U.S.A.