

Clean & Simple™ - Ultrasonic Cleaning Systems

Tuttnauer™
Your Sterilization & Infection Control Partners

Clean &
Simple



- High efficiency 60 kHz of cavitation for a thorough and quieter cleaning
- Durable vinyl clad steel housing... easy to clean and resists scratching
- Stainless steel tank with cover
- Full line of accessories available



Clean & Simple™ - Ultrasonic Cleaning Systems

Tuttnauers' Clean & Simple Ultrasonic Cleaning Systems are manufactured with only the highest quality components. These machines provide 60 kHz of quiet cleaning power which is less aggressive to instruments while still removing all surface contamination. The contemporary design and quiet operation compliment any office environment. Clean & Simple Ultrasonic Systems are available in 2 sizes with a full selection of accessories to meet your specific needs. The 1 gallon holds a standard selection of instruments or 2 half size cassettes. The 3 gallon holds 3 full size cassettes.



Available accessories

System Specifications:		CSU1	CSU3
Tank Dimensions Inches (Cms.)	W	9 3/8 (23.8)	15 1/2 (39.4)
	H	4 (10.2)	6 (15.2)
	D	5 3/8 (13.7)	9 (22.9)
External Dimensions Inches (Cms.)	W	12 5/8 (32)	19 (48.3)
	H	8 3/4 (22.2)	13 1/2 (34.3)
	D	9 7/8 (25)	16 (40.6)
Tank Capacity Gallons (Liters)		0.8 (3.4)	3 (12.7)



Tuttnauer Clean & Simple™ Ultrasonic Tablets

For ultrasonic cleaning, instrument pre-soak and evacuation systems

- Two tablets added to warm water makes 1 gallon ultrasonic / enzymatic solution
- Breaks down protein and organic matter such as blood, mucous and tissue
- Anti-corrosive to protect instruments
- Can be used with any ultrasonic cleaning system • Cleans and deodorizes all evacuation equipment • Non-toxic, bio-degradable, non-caustic • Eliminates heavy handling, mixing and storage of liquid concentrates • Available in 2 sizes: Box of 64 tablets or Box of 144 tablets.



Fully Automatic
EZ Series



Manual
M Series



Valueklave 1730

Tuttnauer Tabletop Autoclaves

Tuttnauer USA Co., Ltd., 25 Power Drive, Hauppauge, NY 11788
Tel: (800) 624 5836, (631) 737 4850, Fax: (631) 737 0720
Email: info@tuttnauerUSA.com, www.tuttnauerUSA.com



LTR330-0009 / ULF / USA / 05.07