



Ultra Procedure Chair with Wheelbase®, Programmable Hand and Foot Controls: Model 4011-650-300

CHAIR SPECIFICATIONS

- Height:** Min. 18.5" (± .5") 46.9cm (±1.27cm) at seat
Max. 37" at seat (93.9cm)
- Seat Width:** 28" (71.12cm)
- Back Width:** 24.5" (62.23cm)
- Length:** 57" – 76" with foot section extended
(185.42cm – 205.74cm)
- Leg rest:** 28" W x 13" D (71.12cm W x 33.02cm D)
- OneTouch WheelBase®:** Integrated Central-Locking
Caster System
- Stirrups:** Four Lateral Positions, Adjustable Length
- Dual Electro-Mechanical Actuator System** for Efficient
Speed
- Electrical:** 120 VAC, 60HZ
- International Model 4011-650-301:** 220/240V, 50/60Hz

ACCESSORIES

- 044 - Adjustable/Self Leveling Arms Pair- Field Installable
- 050 - Removable U-shaped Headrest
- 051 - Removable Rectangular Headrest
- 251 - Bierhoff Knee Crutches
- 467 - Mounting Plate for Welch Allyn Light Bracket #48955 Patient Right
- 468 - Mounting Plate for Welch Allyn Light Bracket #48955 Patient Left
- 830 - Patient Safety Straps
- 900 - Lifting Bracket Kit
- 920 - Easy Installation Kit
- Fire retardant foam (CAL133) available

COLORS



KEY FEATURES

- Chair Designed for Universal Use
- OneTouch Wheelbase®
- 650 lb. (294.84kg) Patient Weight Capacity
- Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
- Hand Control with OneTouch Patient Positioning®
2 Programmable Positions & Return-to-Home
- Foot Control with OneTouch Patient Positioning®
2 Programmable Positions & Return-to-Home
- Removable Seamless Upholstered Cushion
- Double-Articulating Headrest with Adjustable Height
- Interchangeable Headrest Options
- Stirrups
- Stainless Steel Drain Pan
- Trendelenburg, Lithotomy, Semi-Fowler Positions
- Hospital-Grade Electrical Receptacle with Child Safety Cover
(Patient Left)
- Concealed Paper Roll Holder
- 3 Year Limited Warranty

