

Two Function Bariatric Power Table



Standard Features

- Ultra comfort Seamless upholstery with antimicrobial guard and built-in pillow
- ADA Compliant with 19" low access
- All-steel heavy duty construction
- Articulating power backrest
- Power hi-lo
- Two function hand control
- Concealed paper roll holder (18" – 21") (45.7 cm – 53.3cm)
- Large, full extension storage drawer on roller bearing drawer guides
- 10" (25.4cm) leg section
- Electro-mechanical actuator system
- 800 lb. (363 kg) weight capacity
- 3 year limited warranty

Optional Accessories

- **035** - Armboard assembly (specify patient right or left/factory installed)
- **046** - Accessory rails (pair/specify head or seat section/factory installed)
- **450** - Two function foot control (substitution for hand control/64-)
- **451** - Simultaneous two function foot and hand control (69-)
- **802** - Paper cutter band
- **810** - Restraint straps assembly
- **856**- 220/240 Volt (factory installed)
- **5005F** – Fire retardant foam

Table Specifications

- **Height:**
Min. 19" (48.3cm) at seat
Max. 35" (88.9cm) at seat
- **Width:** 30" (76.2cm)
- **Length:** 70" (177.8cm) to 80" (203.2cm) with leg section extended
- **Storage drawer:** 16" D x 19" W x 4.5" H (40.6cm D x 48.3cm W x 11.4cm H)
- **Electrical:** 120 VAC, 60 HZ or 220/240VAC, 60HZ

Shown with -046 accessory rail

Color Chart and Part Number:

5005RI River Rock	5005SN Sand Grey	5005MO Mocha	5005AD Adobe
5005ST Steel Blue	5005MI Midnight Blue	5005SG Sage	5005ON Onyx