

Manual Exam Table with Four Storage Drawers



Color Chart and Part Number:

5240RI River Rock	5240SN Sand Grey	5240MO Mocha	5240AD Adobe
5240ST Steel Blue	5240MI Midnight Blue	5240SG Sage	5240ON Onyx

Standard Features

- 500 Pound capacity (226.8kg)
- Removable ultra comfort seamless top
- 5 year limited warranty
- Galvanized steel base, 18 gauge, powder coated
- Easy-reach articulating back adjustment
- One touch, easy-glide stirrups with 3 lateral positions and friction lock length adjustment
- Slide out leg section with removable upholstered leg rest
- Removable 9"x9" plastic drip pan (22.9cmx22.9cm)
- Pull-out all steel non slip step 7" high (17.8cm)
- 21" Paper roll holder (53.3 cm)
- 2 front drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cmx39.4cmx10.2cm)
- 2 reversible side storage drawers on full extension roller bearing glides (18.25"x15.5"x4" each or 46.4cmx39.4cmx10.2 cm)

Optional Accessories

- **010** - Leg extension booster pad
- **037** - Armboard assembly (specify patient right or left side / factory installed)
- **045** - Accessory rail (pair / specify head or seat section / factory installed)
- **251** - Bierhoff knee crutches (pair)
- **350** - Pelvic tilt
- **365** - I.V. pole (factory installed)
- **375** - Safety rails assembly (pair / factory installed)
- **701** - Removable standard premium top
- **701F** - Removable standard premium top with CAL-133 foam
- **810** - Restraint straps assembly (factory installed)
- **820** - Stainless steel drain pan (factory installed)
- **822** - Stainless steel drain pan with spout (factory installed)
- **849** - Drawer warmer 220V (factory installed)
- **851** - Drawer warmer 110V (factory installed)
- **852** - Hospital grade receptacle 110V with child safe cover (factory installed / patient left side)
- **855** - Hospital grade receptacle 220V with child safe cover (factory installed / patient left side)
- **5240F** - Fire retardant foam
- **5240B** - Base only
- SilverSan antimicrobial powder coated base
- Welch Allyn light bracket

Table Specifications

- **Height:** 31.8" (80.8cm)
- **Length:** 57.5" - 74" (146.1cm-188cm)
- **Width:** 27" (68.6cm)
- **Stirrups:** Infinitely adjustable