

2 Drawers & 2 Doors, 24" Base Cabinet



- All-Steel Construction
- 2 Heavy Duty All-Steel Doors
- Commercial Grade Hinges
- 2 Spacious All-Steel Drawers
- Heavy Duty Roller Bearing Drawer Guides
- Countertop Included
- Available In All 12 UMF Medical Colors

The UMF Medical 6014 Modular Base Cabinet features 24" (61.0cm) high pressure laminate countertop, two drawers and two doors.

3 year warranty included.

Color Chart:



Standard Features

- 20 gauge steel construction
- High pressure laminate countertop with 4" (10.2cm) backsplash
- Full extension roller bearing drawer glides
- Commercial grade, adjustable, spring assisted door hinges
- Interchangeable drawer and door fronts
- 12 accent colors available
- 3 year limited warranty

Specifications

Base cabinet only (including drawers):

- Width: 24" (61.0cm)
- Height: 34.5" (87.6cm)
- Depth: 18" (45.7cm)

Base cabinet with countertop:

- Width: 24" (61.0cm)
- Height: 35.875" (91.1cm) – does not include height of backsplash
- Height: 39.75" (101.0cm) – includes height of backsplash
- Depth: 19" (48.3cm)

Options

- 502 - Double door lock
- 505 - Single drawer lock
- 507 - Additional keys (x2)