

Large Instrument/ Storage and Supply Cabinet

Model 7142

The UMF Medical 7142 Large Instrument / Storage and Supply Cabinet has two lockable upper doors with glass panes, double drawers and double doors in the lower section, three adjustable glass shelves in the upper section, and one removable, non-adjustable metal shelf in the lower section. Comes with a 3-year limited warranty. Available in all UMF Medical accent colors.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › 20 gauge steel construction
- › Standard lock on upper right door
- › Adjustable leveling glides
- › Removable, non-adjustable shelf in lower section
- › Safety release lever on upper left door
- › Door and drawer fronts available in all accent colors
- › Heavy-duty concealed door hinges
- › 3 year limited warranty

Specifications

Width: 31.5" (80cm)
 Height: 66.5" (168.9cm)
 Depth: 15.5" (39.4cm)
 Net Weight: 172lb (78kg)
 Shipping Weight: 248lb (112.5kg)
 Shipping Dimensions: L48xW40xH77
 (121.9cm x 101.6cm x 195.6cm)

Accessories

070 – 110 Volt Light
 071 – 220 Volt Light
 502 – Double door lock (factory installed)
 505 – Single drawer lock (factory installed)
 507 – Additional keys (x2)
 025 – Custom door/drawer front accent color